碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoIHC™ S100 Rabbit Monoclonal Antibody

产品编号	产品名称	包装
AG8670	BeyoIHC™ S100 Rabbit Monoclonal Antibody	50µl

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IHC-P, Immunohistochemistry-paraffin; IHC-Fr, Immunohistochemistry-frozen; IF, Immunofluorescence; ICC, Immunocytochemistry; WB, Western blot; FC, Flow Cytometry.

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About this Antibody	
Name	BeyoIHC™ S100 Rabbit Monoclonal Antibody
Category	Recombinant Rabbit Monoclonal Antibody (mAb); Primary antibody; IHC抗体
Isotype	IgG
Clone number	BA1263
Purification method	Protein A affinity purification
Positive samples	-
Cellular location	Cytoplasm, Nuclears
Customer validation	-
Epitope retrieval	Heat induced epitope retrieval with Tris-EDTA buffer (pH 9.0).
About the Immunogen	
Immunogen	Synthetic peptide corresponding to residues within aa50-150 of S100 was used as an immunogen.
Sequence	-
Gene ID	6285

Swiss Prot	P04271
Synonyms	NEF; S100; S100-B; S100beta; S100B
Category	Cell Cycle
Background	S100 is a 21kDa highly acidic and water soluble protein, which is expressed in a wide variety of normal and neoplastic cells of mesodermal, neuroectodermal, and epithelial origin. The protein is a dimer, having two subunits, alpha and beta, with extensive sequence homology. There are 3 forms of S-100 protein: alpha-alpha, known as S-100A0; alpha-beta, known as S-100A and beta-beta, known as S-100B. The immunohistochemical evaluation of S-100 (beta) protein expression is important in the diagnosis of undifferentiated malignant tumours of unknown primary origin and should be included in the so-called primary panel. S100 is a very sensitive marker for malignant melanoma of all types. S100 may be used in the differential diagnosis of sarcomas and spindle cell tumours.

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c. 一抗孵育(Primary antibody incubation)

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f. 镜检拍照(Microscopic examination and photography)

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
图1. Immunohistochemistry of paraffin-embedded human chordoma using S100 antibody at dilution of 1:200 (200x lens).

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1. Weber J, Peng H, Rader C. Exp Mol Med. 2017. 49(3):e305.
2. Seeber S, Ros F, Thorey I, Tiefenthaler G, Kaluza K, et al. PLoS One. 2014. 9(2):e86184.
3. Zhang Z, Liu H, Guan Q, Wang L, Yuan H. Front Immunol. 2017. 8:494.

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