碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

Poly-L-lysine溶液

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1. Song CY, Wang ZY, Yang J, Zhang RH, Cui YP. Effects of solid substrate on the SERS - based immunoassay: a comparative study. Journal of Raman Spectroscopy. 2011 March;42(3):313-8.
2. Liang GJ, Zhang XF, Wang JJ, Sun YC, Sun XF, Cheng SF, Li L, De Felici M, Shen W. Activin A Accelerates the Progression of Fetal Oocytes Throughout Meiosis and Early Oogenesis in the Mouse. Stem Cells Dev. 2015 Oct 15;24(20):2455-65.
3. Shen X, Ma L, Dong W, Wu Q, Gao Y, Luo C, Zhang M, Chen X, Tao L. Autophagy regulates intracerebral hemorrhage induced neural damage via apoptosis and NF-κB pathway. Neurochem Int. 2016 Jun;96:100-12.
4. Mu Zhang, Fei Qian, Qingling Liu, Cheng Qian, Pyone Myat Thu, Yanyan Wang, Zu-Guo Zheng, Hua Yang, Ping Li and Xiaojun Xu. Evaluation of structure-activity relationships of ginsenosides against amyloid β induced pathological behaviours in transgenic Caenorhabditis elegans. Rsc Adv 2017; 7: 40095-40104..
5. Cheng W, Yu P, Wang L, Shen C, Song X, Chen J, Tang F, Yang Q. Sonic hedgehog signaling mediates resveratrol to increase proliferation of neural stem cells after oxygen-glucose deprivation/reoxygenation injury in vitro. Cell Physiol Biochem. 2015;35(5):2019-32.
6. Shen X, Ma L, Dong W, Wu Q, Gao Y, Luo C, Zhang M, Chen X, Tao L. Autophagy regulates intracerebral hemorrhage induced neural damage via apoptosis and NF-κB pathway. Neurochem Int. 2016 Jun;96:100-12
7. Zhang J, Li M, Li Z, Shi J, Zhang Y, Deng X, Liu L, Wang Z, Qi Y, Zhang H. Deletion of the Type IV Secretion System Effector VceA Promotes Autophagy and Inhibits Apoptosis in Brucella-Infected Human Trophoblast Cells. Curr Microbiol. 2019 Apr;76(4):510-519.