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1. Wang Y, Deng L, Zhong H, Wang Y, Jiang X, Chen J. Natural plant extract tubeimoside I promotes apoptosis-mediated cell death in cultured humanhepatoma (HepG2) cells. *Biol Pharm Bull.* 2011;34(6):831-8.
2. Yanchun Wei , Qun Chen , Baoyan Wu , Aiguo Zhou and Da Xing. High-sensitivity in vivo imaging for tumors using a spectral up-conversion nanoparticle NaYF<sub>4</sub>: Yb<sup>3+</sup>, Er<sup>3+</sup> in cooperation with a microtubulin inhibitor. *Nanoscale.* 2012 May;13:3901-9.
3. Zheng M,Zheng Y,Xie L, Chang W, Gu N Ji M. Orally Active Aurora A/B Kinase Inhibitor, AM - 005, Suppresses the Growth of Human Colon Carcinoma Cells. *Drug Development Research.* 2013 Jun;74(4):272-81.
4. Feng LX, Li M, Liu YJ, Yang SM, Zhang N. Synergistic enhancement of cancer therapy using a combination of ceramide and docetaxel. *Int J Mol Sci.* 2014 Mar 10;15(3):4201-20.
5. Jin H, Liang Q, Chen T, Wang X. Resveratrol protects chondrocytes from apoptosis via altering the ultrastructural and biomechanical properties: an AFM study.*PLoS One.* 2014 Mar 14;9(3):e91611.
6. Lu Z, Li H, Hou C, Peng Y, Long J, Liu J. Endogenously generated amyloid- $\beta$  increases stiffness in human neuroblastoma cells. *Eur Biophys J.* 2016 Nov 16. [Epub ahead of print]
7. Tu LY, Bai HH, Cai JY, Deng SP. The mechanism of kaempferol induced apoptosis and inhibited proliferation in human cervical cancerSiHa cell: From macro to nano. *Scanning.* 2016 Nov;38(6):644-653.
8. Liu L, Hu Y, Shen YF, Wang GX, Zhu B. Evaluation on antiviral activity of coumarin derivatives against spring viraemia of carp virus in epithelioma papulosum cyprini cells. *ANTIVIR RES.* 2017

- Aug;144:173-185.
9. Zhang C,Wang W,He W,Xi N,Wang Y,Liu L. Dynamic Model for Characterizing Contractile Behaviors and Mechanical Properties of a Cardiomyocyte. *Biophys J* . 2018 Jan 9;114(1):188-200.
  10. Xue YZB,Niu YM,Tang B,Wang CM. PCL/EUG scaffolds with tunable stiffness can regulate macrophage secretion behavior. *PROG BIOPHYS MOL BIO* . 2018 May 21. pii: S0079-6107(17)30320-6.
  11. Zhao YQ,Mu DL,Wang D,Han YL,Hou CC,Zhu JQ. Analysis of the function of KIF3A and KIF3B in the spermatogenesis in *Boleophthalmus pectinirostris*. *Fish Physiol Biochem* . 2018 Jun;44(3):769-788.
  12. Zhang MQ,Li JR,Peng ZG,Zhang JP. Differential Effects of Autophagy-Related 10 Protein on HCV Replication and Autophagy Flux Are Mediated by Its Cysteine44 and Cysteine135. *Front Immunol* . 2018 Sep 24;9:2176.
  13. Shen YF,Liu L,Feng CZ,Hu Y,Chen C,Wang GX,Zhu B. Synthesis and antiviral activity of a new coumarin derivative against spring viraemia of carp virus. *FISH SHELLFISH IMMUN* . 2018 Oct;81:57-66.
  14. Liu L,Zhang C,Wang W,Xi N,Wang Y. Regulation of C2C12 Differentiation and Control of the Beating Dynamics of Contractile Cells for a Muscle-Driven Biosyncretic Crawler by Electrical Stimulation. *Soft Robot* . 2018 Oct 18.
  15. Lu X,Wang Q,Gu H,Zhang X,Qi Y,Liu Y. Whole exome sequencing identified a second pathogenic variant in HOMER2 for autosomal dominant non-syndromic deafness. *Clin Genet* . 2018 Nov;94(5):419-428.
  16. Chen WC,Liu L,Shen YF,Hu Y,Ling F,Wang GX,Zhu B. A new coumarin derivative plays a role in rhabdoviral clearance by interfering glycoprotein function during the early stage of viral infection. *Cell Signal* . 2018 Nov;51:199-210.
  17. Li Z,Liu X,Ma J,Zhang T,Gao X,Liu L. hnRNPK modulates selective quality-control autophagy by downregulating the expression of HDAC6 in 293 cells. *Int J Oncol* . 2018 Nov;53(5):2200-2212.
  18. Gao Y,Elamin E,Zhou R,Yan H,Liu S,Hu S,Dong J,Wei M,Sun L,Zhao Y. FKBP51 promotes migration and invasion of papillary thyroid carcinoma through NF- $\kappa$  B-dependent epithelial-to-mesenchymal transition. *Oncol Lett* . 2018 Dec;16(6):7020-7028.
  19. Yilu Su,Zhi Zhou,Xiaopeng Yu. Possible roles of glutamine synthetase in responding to environmental changes in a scleractinian coral. *Mol Biol Rep*. 2018 Dec;45(6):2115-2124.;doi: 10.1007/s11033-018-4369-3.
  20. Luo Z,Xu X,Sho T,Zhang J,Xu W,Yao J,Xu J. ROS-induced autophagy regulates porcine trophectoderm cell apoptosis, proliferation, and differentiation. *AM J PHYSIOL-CELL PH* . 2019 Feb 1;316(2):C198-C209.
  21. Tang JC,Wu K,Zheng X,Xu M,Dai Y,Wei SS,Cai XJ. GSK923295 as a potential antihepatocellular carcinoma agent causing delay on liver regeneration after partial hepatectomy. *CHINESE MED J-PEKING*. 2019 Feb 5;132(3):311-318.
  22. Sun J,Zhou Q,Hu X. Integrating multi-omics and regular analyses identifies the molecular responses of zebrafish brains to graphene oxide: Perspectives in environmental criteria. *ECOTOX ENVIRON SAFE*. 2019 Sep 30;180:269-279.
  23. Wei YL,Yang T,Kovacs T,Yang WX. C-terminal kinesin motor es-KIF1C1 regulates nuclear formation during spermiogenesis in Chinese mitten crab *Eriocheir sinensis*. *Gene*. 2019 Nov 30;719:144074.
  24. Xiaolei Wang,Kaiyan Lou,Xiaodong Song,Huijuan Ma,Xinyi Zhou,Huan Xu,Wei Wang. Mebendazole is a potent inhibitor to chemoresistant T cell acute lymphoblastic leukemia cells. *TOXICOL APPL PHARM* . 2020 Jun 1;396:115001.;doi: 10.1016/j.taap.2020.115001.
  25. Yiming Niu,Lintao Wang,Na Yu,Panfei Xing,Zhenzhen Wang,Zhangfeng Zhong,Yanxian Feng,Lei Dong,Chunming Wang. An "all-in-one" scaffold targeting macrophages to direct endogenous bone repair in situ. *Acta Biomater*. 2020 Jul 15;111:153-169.;doi: 10.1016/j.actbio.2020.05.023.
  26. Zhen Chen,Chao Chen,Tingting Zhou,Chao Duan,Qianqian Wang,Xiaohui Zhou,Xia Zhang,Fangrong Wu,Yunfen Hua,Fan Lin. A high-throughput drug combination screen identifies an anti-glioma synergism between TH588 and PI3K inhibitors. *Cancer Cell Int*. 2020 Jul 23;20:337.;doi: 10.1186/s12935-020-01427-0.
  27. Qiangqiang Wang,Jialei Bai,Baoan Ning,Longxing Fan,Tieqiang Xu,Zhiqiang Chen,Shutao Guo. Acid-sensitive PEGylated paclitaxel produg nanoparticles for cancer therapy: Effect of PEG length on antitumor efficacy. *J Control Release*. 2020 Oct 10;326:265-275.;doi: 10.1016/j.jconrel.2020.07.022.
  29. Shuang Li,Zhiqi Wang,Huili Tong,Shufeng Li,Yunqin Yan. TCP11L2 promotes bovine skeletal muscle-derived satellite cell migration and differentiation via FMNL2. *J Cell Physiol*. 2020 Oct;235(10):7183-7193.;doi: 10.1002/jcp.29617.
  30. Yiwen Jiang,Dan Hong,Zhefeng Lou,Xuezi Tu,Longjin Jin. Lupeol inhibits migration and invasion of colorectal cancer cells by suppressing RhoA-ROCK1 signaling pathway. *N-S ARCH PHARMACOL*. 2020 Nov;393(11):2185-2196.;doi: 10.1007/s00210-020-01815-3.
  31. Mengyao Liu,Manfei Fu,Xiaoye Yang,Guoyong Jia,Xiaoqun Shi,Jianbo Ji,Xianghong Liu,Guangxi Zhai. Paclitaxel and quercetin co-loaded functional mesoporous silica nanoparticles overcoming multidrug resistance in breast cancer. *COLLOID SURFACE B*. 2020 Dec;196:111284.;doi: 10.1016/j.colsurfb.2020.111284.
  32. Jiaying Li,Jia Qi,Guangxin Yao,Qinling Zhu,Xinyu Li,Rui Xu,Zhenyi Zhu,Hanting Zhao,Yuan Wang,Ying Ding,Yun Sun. Deficiency of Sirtuin 1 Impedes Endometrial Decidualization in Recurrent Implantation Failure Patients. *Front Cell Dev Biol*. 2021 Jan 28;9:598364.;doi: 10.3389/fcell.2021.598364.
  33. Hao Cai,Yufan Xiang,Yujun Zeng,Zhiqian Li,Xiuli Zheng,Qiang Luo,Hongyan Zhu,Qiyong Gong,Zhongwei Gu,Yanhui Liu,Hu Zhang,Kui Luo. Cathepsin B-responsive and gadolinium-labeled branched glycopolymer-PTX conjugate-derived nanotheranostics for cancer treatment. *Acta Pharm Sin B*. 2021 Feb;11(2):544-559.;doi: 10.1016/j.apsb.2020.07.023.

Version 2021.07.02