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- Lin Zhang, Mei Xu, Qiao Ren, Gang Liu, Shulin Meng, Kang Xiahou, Yongxiang Zhang, Ning Jiang, Wenxia Zhou. Human Induced Pluripotent Stem Cell-Derived Neural Cells from Alzheimer's Disease Patients Exhibited Different Susceptibility to Oxidative Stress. *Stem Cells Dev.* 2020 Nov;15;29(22):1444-1456.
- Duo Wang, Jun Zhou, Weimin Fang, Cuiqing Huang, Zerong Chen, Meng Fan, Ming-Rong Zhang, Zeyu Xiao, Kuan Hu, Liangping Luo. A multifunctional nanotheranostic agent potentiates erlotinib to EGFR wild-type non-small cell lung cancer. *Bioact Mater.* 2021 Nov 4:13:312-323.
- Yang Zhang, Shan Liu, Jiu-Long Ma, Chen Chen, Peng Huang, Jia-Hua Ji, Di Wu, Li-Qun Ren. Apocynum venetum leaf extract alleviated doxorubicin-induced cardiotoxicity through the AKT/Bcl-2 signaling pathway. *Phytomedicine.* 2022 Jan;94:153815.
- Haiwei Ma, Chenglong Xie, Zhengtai Chen, Gaolu He, Zihan Dai, Hanchen Cai, Haojie Zhang, Hongwei Lu, Hongqiang Wu, Xinli Hu, Kailiang Zhou, Gang Zheng, Huazi Xu, Cong Xu. MFG-E8 alleviates intervertebral disc degeneration by suppressing pyroptosis and extracellular matrix degradation in nucleus pulposus cells via Nrf2/TXNIP/NLRP3 axis. *Cell Death Discov.* 2022 Apr 19;8(1):209.
- Xuehan Sang, Xiuhong Zhao, Lianqi Yan, Xing Jin, Xin Wang, Jianjian Wang, Zhenglu Yin, Yuxin Zhang, Zhaoxiang Meng. Thermosensitive Hydrogel Loaded with Primary Chondrocyte-Derived Exosomes Promotes Cartilage Repair by Regulating Macrophage Polarization in Osteoarthritis. *Tissue Eng Regen Med.* 2022 Jun;19(3):629-642.
- Xu Peng, Can Cheng, Lunli Yue, Yan Liu, Xixun Yu. A Comparative Study Between Porcine Peritoneum and Pericardium as Cardiovascular Material. *Tissue Eng Part C Methods.* 2022 Jun;28(6):272-284.
- Yue Zhang, Ying Wang, Lin Zeng, Yixiang Liu, Huijing Sun, Shanshan Li, Siyu Wang, Longjun Shu, Naixin Liu, Saige Yin, Junsong Wang, Dan Ni, Yutong Wu, Ying Yang, Li He, Buliang Meng, Xinwang Yang.

- Amphibian-derived peptide homodimer OA-GL17d promotes skin wound regeneration through the miR-663a/TGF- β 1/Smad axis. Burns Trauma. 2022 Jul 12:10:tkac032.
8. Yongdie Deng, Fenghua Zhang, Menglu Jiang, Yanju Liu, Huiping Yuan, Jinsong Leng. Programmable 4D Printing of Photoactive Shape Memory Composite Structures. ACS Appl Mater Interfaces. 2022 Sep 21;14(37):42568-42577.
 9. Yi-Dan Zheng, Li-Mei Ma, Jin-Jian Lu, Tian Chai, Mohammad Reza Naghavi, Jun-Yi Ma, Chun-Yan Sang, Jun-Li Yang. Nardoguajanone L Isolated from Nardostachys jatamansi Improved the Effect of Gemcitabine Chemotherapy via Regulating AGE Signaling Pathway in SW1990 Cells. Molecules. 2022 Oct 13;27(20):6849.
 10. Changgui Huang, Jielin Zhou, Jin Rao, Xinyi Zhao, Xiumei Tian, Fupo He, Haishan Shi. Fabrication of strontium carbonate-based composite bioceramics as potential bone regenerative biomaterials. Colloids Surf B Biointerfaces. 2022 Oct:218:112755.
 11. Cunqing Kong, Banghao Xu, Guanhua Qiu, Meng Wei, Mengqi Zhang, Shengxian Bao, Jiali Tang, Lequn Li, JunJie Liu. Multifunctional Nanoparticles-Mediated PTT/PDT Synergistic Immune Activation and Antitumor Activity Combined with Anti-PD-L1 Immunotherapy for Breast Cancer Treatment. Int J Nanomedicine. 2022 Nov 17:17:5391-5411.
 12. Runjie Guo, Ping Chen, Tiantian Fu, Ren Zhang, Yuan Zhu, Nange Jin, Hong Xu, Yong Xia, Xuesong Tian. Activation of Delta-Opioid Receptor Protects ARPE19 Cells against Oxygen-Glucose Deprivation/Reoxygenation-Induced Necroptosis and Apoptosis by Inhibiting the Release of TNF- α . J Ophthalmol. 2022 Nov 19:2022:2285663.
 13. Yan Li, Hong-Ying Li, Jun Shao, Lingpeng Zhu, Tian-Hua Xie, Jiping Cai, Wenjuan Wang, Meng-Xia Cai, Zi-Li Wang, Yong Yao, Ting-Ting Wei. GRP75 Modulates Endoplasmic Reticulum-Mitochondria Coupling and Accelerates Ca²⁺-Dependent Endothelial Cell Apoptosis in Diabetic Retinopathy. Biomolecules. 2022 Nov 29;12(12):1778.
 14. Nana Qiao, Jinghui Zhang, Ya Zhang, Xinjie Liu. Synergistic regulation of microglia differentiation by CD93 and integrin β 1 in the rat pneumococcal meningitis model. Immunol Lett. 2022 Dec:251-252:63-74.
 15. Yiting Lei, Yuping Wang, Jieliang Shen, Zhengwei Cai, Chen Zhao, Hong Chen, Xiaoji Luo, Ning Hu, Wengu Cui, Wei Huang. Injectable hydrogel microspheres with self-renewable hydration layers alleviate osteoarthritis. Sci Adv. 2022 Feb 4;8(5):eabl6449.
 16. Chen Zhang, Han Wang, Xinhui Yang, Zi Fu, Xiuru Ji, Yifan Shi, Jie Zhong, Weiguo Hu, Youqiong Ye, Zhengting Wang, Dalong Ni. Oral zero-valent-molybdenum nanodots for inflammatory bowel disease therapy. Sci Adv. 2022 Sep 16;8(37):eabp9882.
 17. Huaying He, Weikang Zhou, Jing Gao, Fan Wang, Shaobing Wang, Yan Fang, Yang Gao, Wei Chen, Wen Zhang, Yunxiang Weng, Zhengchao Wang, Haiqing Liu. Efficient, biosafe and tissue adhesive hemostatic cotton gauze with controlled balance of hydrophilicity and hydrophobicity. Nat Commun. 2022 Jan 27;13(1):552

Version 2024.03.12