碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Wen-guang Zhou, Yu-tuo Wei, Kun Huang and Ri-bo Huang. Cloning and expression of glycerol dehydratase gene from Klebsiella pneumoniae in Escherichia coli. Indust Microbiol. 2004;34(2): 1-4.
2. Wang B, Zhang J, Duan C, Wang D. Primary Study on Sound Effecting the Gene Expressions of Arabidopsis thaliana. J Chongqing Univ (Nat Sci Ed). 2005; 8(11): 138-141.
3. Li K, Chen L, Zhou Z, Li E, Zhao X, Guo H. The site of vitellogenin synthesis in Chinese mitten-handed crab Eriocheir sinensis. Comp Biochem Physiol B Biochem Mol Biol. 2006 Apr;143(4):453-8.
4. Yuan DD, Han YH, Wang CG, Chi CW. From the identification of gene organization of alpha conotoxins to the cloning of novel toxins. Toxicon. 2007 Jun 15;49(8):1135-49.
5. Zhou B, Huang C, Yang J, Lu J, Dong Q, Sun LZ. Preparation of heteroduplex enhanced green fluorescent protein plasmid for in vivo mismatch repair activity assay. Anal Biochem. 2009 May 1;388(1):167-9.
6. Yongjin J. Zhou, Fan Yang, Sufang Zhang, Haidong Tan and Zongbao K. Zhao. Efficient gene disruption in Saccharomyces cerevisiae using marker cassettes with long homologous arms prepared by the restriction-free cloning strategy. WORLD JOURNAL OF MICROBIOLOGY AND BIOTECHNOLOGY,2011 Dec; 27(12):2999-3000.
7. Zhou Y, Wang L, Yang F, Lin X, Zhang S, Zhao ZK. Determining the extremes of the cellular NAD(H) level by using an Escherichia coli NAD(+)-auxotrophic mutant. Appl Environ Microbiol. 2011 Sep; 77(17):6133-40.
8. Guo N, Chen R, Li Z, Liu Y, Cheng D, Zhou Q, Zhou J, Lin Q. Hepatitis C virus core upregulates the methylation status of the RASSF1A promoter through regulation of SMYD3 in hilar cholangiocarcinoma cells. Acta Biochim Biophys Sin (Shanghai). 2011 May;43(5):354-61.
9. Zhang CW, Zhang XL, Xia YJ, Cao YX, Wang WJ, Xu P, Che YN, Wu XK, Yi L, Gao Q, Wang Y. Association between polymorphisms of the CYP11A1 gene and polycystic ovary syndrome in Chinese women. Mol Biol Rep. 2012 Aug;39(8):8379-85.
10. Qian J, Wu J, Yao B, Lu Y. Preparation of a polyclonal antibody against hypericin synthase and localization of the enzyme in red-pigmented Hypericum perforatum L. plantlets. Acta Biochim Pol. 2012;59(4):639-45.
11. Zhou H, Chen S, Wang W, Wang Z, Wu X, Zhang Z. Nanog inhibits lipopolysaccharide-induced expression of pro-inflammatory cytokines by blocking NF-κB transcriptional activity in rat primary microglial cells. Mol Med Rep. 2012 Mar;5(3):842-6.
12. Hu D, Ju X, Li L, Hu C, Yan L, Wu T, Fu J, Qin M. Improved in situ saccharification of cellulose pretreated by dimethyl sulfoxide/ionic liquid using cellulase from a newly isolated Paenibacillus sp. LLZ1. Bioresour Technol. 2015 Nov 24;201:8-14.
13. Fu S, Fan L, Pan X, Sun Y, Zhao H. Integrin αv promotes proliferation by activating ERK 1/2 in the human lung cancer cell line A549. Mol Med Rep. 2015 Feb;11(2):1266-71.
14. Hu D, Ju X, Li L, Hu C, Yan L, Wu T, Fu J, Qin M. Improved in situ saccharification of cellulose pretreated by dimethyl sulfoxide/ionic liquid using cellulase from a newly isolated Paenibacillus sp. LLZ1. Bioresour Technol. 2016 Feb;201:8-14.
15. He WT, Zhang LM, Li C, Li SY, Ding ZC, Fang ZM, Meng FY, Chen ZK, Zhou P. Short-term MyD88 inhibition ameliorates cardiac graft rejection and promotes donor-specific hyporesponsiveness of skin grafts in mice. Transpl Int. 2016 Aug;29(8):941-52.
16. Zhou ZJ, Sun L. Edwardsiella tarda-Induced Inhibition of Apoptosis: A Strategy for Intracellular Survival. Front Cell Infect Microbiol. 2016 Jul 14;6:76.
17. Liu X, Nie S, Yu Q, Wang X, Huang D, Xie M. Downregulation of steroid hormone receptor expression and activation of cell signal transduction pathways induced by a chiral nonylphenol isomer in mouse sertoli TM4 cells. Environ Toxicol. 2017 Feb;32(2):469-476.
18. Du X, Shen T, Wang H, Qin X, Xing D, Ye Q, Shi Z, Fang Z, Zhu Y, Yang Y, Peng Z, Zhao C, Lv B, Li X, Liu G, Li X. Adaptations of hepatic lipid metabolism and mitochondria in dairy cows with mild fatty liver. J Dairy Sci. 2018 Oct;101(10):9544-9558.
19. Jin G, Chen Z, Zhang J, Song J, Shi J, Zhou B. Association of brain natriuretic peptide gene polymorphisms with chronic obstructive pulmonary disease complicated with pulmonary hypertension and its mechanism. BIOSCIENCE REP. 2018 Oct 2;38(5). pii: BSR20180905
20. Ma L, Zeng F, Cong F, Huang B, Huang R, Ma J, Guo P. Development of a SYBR green-based real-time RT-PCR assay for rapid detection of the emerging swine acute diarrhea syndrome coronavirus. J Virol Methods. 2019 Mar;265:66-70.
21. Xi Y, Xu CZ, Xie ZZ, Zhu DL, Dong JM. Rapid and visual detection of dengue virus using recombinase polymerase amplification method combined with lateral flow dipstick. MOL CELL PROBE. 2019 Aug;46:101413

Version 2021.09.01