碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

## Annealing Buffer for DNA Oligos(5X)

产品编号	产品名称	包装
D0251	Annealing Buffer for DNA Oligos(5X)	1ml

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1. Qu B, Han X, Tang Y, Shen N. A Novel Vector-Based Method for Exclusive Overexpression of Star-Form MicroRNAs. *PLoS One*. 2012;7(7):e41504.
2. Jiang W, Liang P, Wang B, Fang J, Lang J, Tian G, Jiang J, Zhu TF. Optimized DNA extraction and metagenomic sequencing of airborne microbial communities. *Nat Protoc*. 2015 May;10(5):768-79.
3. Huang G, Gao Q, Zhao Y, Dong Z, Li T, Guan X, Jiang J. A novel siRNA validation system for functional screening of effective RNAi targets in mammalian cells and development of a derivative lentivirus delivery system. *Gene*. 2015 Mar 10;558(2):278-86.
4. Yuanyou Li, Ziyang Yin, Yewei Dong, Shuqi Wang, Óscar Monroig, Douglas R Tocher, Cuihong You. Ppar  $\gamma$  Is Involved in the Transcriptional Regulation of Liver LC-PUFA Biosynthesis by Targeting the  $\Delta 6 \Delta 5$  Fatty Acyl Desaturase Gene in the Marine Teleost *Siganus Canaliculatus* *MAR BIOTECHNOL*. 2019 Feb;21(1):19-29.;doi: 10.1007/s10126-018-9854-0
5. Tian M, Qu J, Li P, Bao Y, Liu J, Ding C, Wang S, Li T, Qi J, Yu S. Identification of novel genes essential for *Brucella abortus* to establish infection by signature-tagged mutagenesis. *Vet Microbiol*. 2019 Mar;230:130-137
6. Fang H, Dong Y, Yue X, Hu J, Jiang S, Xu H, Wang Y, Su M, Zhang J, Zhang Z, Wang N, Chen X. The B-box zinc finger protein MdBBX20 integrates anthocyanin accumulation in response to ultraviolet radiation and low temperature. *Plant Cell Environ*. 2019 Jul;42(7):2090-2104
7. Chen H, Wang J, Cui J, Wang C, Liang S, Liu H, Wen J. Negative regulation of bleomycins biosynthesis by ArsR/SmtB family repressor BlmR in *Streptomyces verticillus*. *APPL MICROBIOL BIOT*. 2019 Aug;103(16):6629-6644
8. Chang Y, Deng Y, Li T, Wang J, Wang T, Tan F, Li X, Tian K. Visual detection of porcine reproductive and respiratory syndrome virus using CRISPR-Cas13a. *Transbound Emerg Dis*. 2020 Mar;67(2):564-571
9. Yaokang Wu, Yanfeng Liu, Xueqin Lv, Jianghua Li, Guocheng Du, Long Liu. CAMERS-B: CRISPR/Cpf1 assisted multiple-genes editing and regulation system for *Bacillus subtilis* *Biotechnol Bioeng*. 2020 Jun;117(6):1817-1825.;doi: 10.1002/bit.27322

Version 2021.09.01