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1. Shi Y, Huang H, Deng X, He X, Yang J, Yang H, Shi L, Mei L, Gao J, Zhou N. Identification and functional characterization of two orphan G-protein-coupled receptors for adipokinetic hormones from silkworm *Bombyx mori*. *J Biol Chem*. 2011 Dec 9;286(49):42390-402.
2. Yang J, Huang H, Yang H, He X, Jiang X, Shi Y, Alatangale D, Shi L, Zhou N. Specific activation of the G protein-coupled receptor BNGR-A21 by the neuropeptide corazonin from the silkworm, *Bombyx mori*, dually couples to the G(q) and G(s) signaling cascades. *J Biol Chem*. 2013 Apr 26;288(17):11662-75.
3. Deng X, Yang H, He X, Liao Y, Zheng C, Zhou Q, Zhu C, Zhang G, Gao J, Zhou N. Activation of *Bombyx mori* neuropeptide G protein-coupled receptor A4 via a G α i-dependent signaling pathway by direct interaction with neuropeptide F from silkworm, *Bombyx mori*. *Insect Biochem Mol Biol*. 2014 Feb;45:77-88.
4. Chen X, Zheng C, Qian J, Sutton SW, Wang Z, Lv J, Liu C, Zhou N. Involvement of β -arrestin-2 and clathrin in agonist-mediated internalization of the human cannabinoid CB2 receptor. *Curr Mol Pharmacol*. 2014;7(1):67-80.
5. Yang W, X Dai, M Liu. Development of an Indirect Competitive Enzyme-linked Immunosorbent Assays Method Based on Immunomagnetic-bead for Analyzing *Listeria monocytogenes* in Food. *Am J Food Technol*. 2016;11(3):100-107.
6. Wang T, Yang Z, Zhou N, Sun L, Lv Z, Wu C. Identification and functional characterisation of 5-HT4 receptor in sea cucumber *Apostichopus japonicus*(Selenka). *Sci Rep*. 2017 Jan 6;7:40247.
7. Shen Z, Chen Y, Hong L, Cui Z, Yang H, He X, Shi Y, Shi L, Han F, Zhou N. BNGR-A25L and -A27 are two functional G protein-coupled receptors for CAPA periviscerokinins in the silkworm *Bombyx mori*. *J Biol Chem*. 2017 Oct 6;292(40):16554-16570.
8. Shen Z, Yang X, Chen Y, Shi L. CAPA periviscerokinins-mediated activation of MAPK/ERK signaling through Gq-PLC-PKC-dependent cascade and reciprocal ERK activation-dependent internalized kinetics of Bom-CAPA-PVK receptor 2. *INSECT BIOCHEM MOLEC*. 2018 Jul;98:1-15.
9. Wang T, Liang J, Xiang X, Chen X, Zhang B, Zhou N, Huang W, Yang J. Pharmacological characterization, cellular localization and expression profile of NPY receptor subtypes Y2 and Y7 in large yellow croaker, *Larimichthys crocea*. *COMP BIOCHEM PHYS B*. 2019 Dec;238:110347.

Version 2021.09.01