碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

BeyoFast™ SYBR Green qPCR Mix (2X, High ROX)

产品编号	产品名称	包装
D7265-1ml	BeyoFast™ SYBR Green qPCR Mix (2X, High ROX)	1ml
D7265-5ml	BeyoFast™ SYBR Green qPCR Mix (2X, High ROX)	5ml
D7265-25ml	BeyoFast™ SYBR Green qPCR Mix (2X, High ROX)	25ml

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1. Zheng GH, Shan Q, Mu JJ, Wang YJ, Zhang ZF, Fan SH, Hu B, Li MQ, Xie J, Chen P, Wu DM, Lu J, Zheng YL. Purple Sweet Potato Color Attenuates Kidney Damage by Blocking VEGFR2/ROS/NLRP3 Signaling in High-Fat Diet-Treated Mice. *Oxid Med Cell Longev*. 2019 Jan 22; 2019:5189819.
2. Cao L, Zhang X, Yuan S, Cheng K, Zhang X. Autophagy induced by enterovirus 71 regulates the production of IL-6 through the p38MAPK and ERK signaling pathways. *MICROB PATHOGENESIS*. 2019 Jun; 131:120-127.
3. Song Meng, Ponnie Robertlee Dolo, Pei Guo, Jian Hong, Chao Li, Xiaocheng Zhu, Dinghua Zhou. The expression of long non-coding RNA LINC01279 in gastric adenocarcinoma and its clinical significance. *Asian J Surg*. 2022 Jun; 45(6):1231-1236.
4. Bailing Zhang, Honglang Li, Jieling Zhang, Yaping Hang, Yi Xu. Activating transcription factor 3 protects alveolar epithelial type II cells from Mycobacterium tuberculosis infection-induced inflammation. *Tuberculosis (Edinb)*. 2022 Jul; 135:102227.
5. Zhen-Yu She, Meng-Fei Xu, Sun-Ying Jiang, Ya-Lan Wei. Kinesin-7 CENP-E is essential for chromosome alignment and spindle assembly of mouse spermatocytes. *Biochim Biophys Acta Mol Cell Res*. 2022 Sep; 1869(9):119306.
6. Haocun Kong, Luxi Yu, Caiming Li, Xiaofeng Ban, Zhengbiao Gu, Zhao Feng Li. Short-Clustered Maltodextrin Activates Ileal Glucose-Sensing and Induces Glucagon-like Peptide 1 Secretion to Ameliorate Glucose Homeostasis in Type 2 Diabetic Mice. *J Agric Food Chem*. 2022 Oct 5; 70(39):12604-12619.
7. Wen-Xin Yu, Yu-Kun Li, Meng-Fei Xu, Chen-Jie Xu, Jie Chen, Ya-Lan Wei, Zhen-Yu She. Kinesin-5 Eg5 is essential for spindle assembly, chromosome stability and organogenesis in development. *Cell Death Discov*. 2022 Dec 13; 8(1):490.
8. Qiang Sun, Qiu Tao, Tianqi Ming, Shun Tang, Hui Zhao, Maolun Liu, Han Yang, Shan Ren, Jiarong Lei, Yuanjing Liang, Yuhui Peng, Minmin Wang, Haibo Xu. Berberine is a suppressor of Hedgehog signaling cascade in colorectal cancer. *Phytomedicine*. 2023 Jun; 114:154792.
9. Xiaofang Zhu, Jing Wan, Xu You, Wanli Yang, Lei Zhao. Circular non-coding RNA circ_0072088 serves as a ceRNA, targeting the miR-1225-5p/WT1 axis to regulate non-small cell lung cancer cell malignant behavior. *Thorac Cancer*. 2023 Jul; 14(20):1969-1979.

Version 2024.03.12