碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Fengxiang Wang, Emily T Graham, Nathchar Naowarajna, Zhennan Shi, Yuqi Wang, Guanglei Xie, Lili Zhou, Wendy Salmon, Jie-Min Jia, Xi Wang, Yuwei Huang, Stuart L Schreiber, Yilong Zou . PALP: A rapid imaging technique for stratifying ferroptosis sensitivity in normal and tumor tissues in situ Cell Chem Biol. 2022 Jan 20;29(1):157-170.e6.
2. Deng Z, Wang Z, Jin J, Wang Y, Bao N, Gao Q, Zhao J . SIRT1 protects osteoblasts against particle-induced inflammatory responses and apoptosis in aseptic prosthesis loosening. Acta Biomater. 2017 Feb;49:541-554.
3. Wang F, Liang YJ, Wu XP, Su XD, Fu LW . ABCG2-overexpressing S1-M1-80 cell xenografts in nude mice keep original biochemistry and cell biological properties. Chin J Cancer. 2012 Mar;31(3):150-8.
4. Deng Z, Jin J, Wang Z, Wang Y, Gao Q, Zhao J . The metal nanoparticle-induced inflammatory response is regulated by SIRT1 through NF- κ B deacetylation in aseptic loosening. INT J NANOMED. 2017 May 10;12:3617-3636.
5. Liu X, Zhang Z, Li P, Zhu L, Wang Y, Wang C . Polypeptide from Chlamys farreri modulates UVB-induced activation of NF-kappaB signaling pathway and protection HaCaT cells from apoptosis. REGUL PEPTIDES. 2009 Feb 25;153(1-3):49-55.
6. Liu HT, Li WM, Li XY, Xu QS, Liu QS, Bai XF, Yu C, Du YG . Chitosan oligosaccharides inhibit the expression of interleukin-6 in lipopolysaccharide-induced human umbilical vein endothelial cells through p38 and ERK1/2 protein kinases. BASIC CLIN PHARMACOL. 2010 May;106(5):362-71.
7. Wang XH, Fu ZG, Zhao Y, Shen W, Wu XY, Wang CY . Survivin: an overexpression protein with notable cellular localization and multiple roles in colon cancer. 《中德临床肿瘤学杂志》. 2010 Sep;9(9):519-23.
8. Wang X, Fu Z, Zhao Y, Wu X, Shen W . Profile of protein expression of the colon cancer cell line SW480 with survivin/shRNA. Eur J Cancer Prev. 2011 May;20(3):190-8.
9. Wang X, Zhao Y, Fu Z, He Y, Xiang D, Zhang L . Prelining autogenic endothelial cells in allogeneic vessels inhibits thrombosis and intimal hyperplasia: an efficacy study in dogs. J Surg Res. 2011 Jul;169(1):148-55.
10. Zhou L, Xue H, Wang Z, Ni J, Yao T, Huang Y, Yu C, Lu L . Angiotensin-(1-7) attenuates high glucose-induced proximal tubular epithelial-to-mesenchymal transition via inhibiting ERK1/2 and p38 phosphorylation. Life Sci. 2012 Mar 10;90(11-12):454-62.
11. Nie J, Cao X, Zhou M, Zhang X, Zhang R, Niu L, Xia Y, Hou L, Wang C . Musca domestica pupae lectin induces apoptosis in HepG2 cells through a NF- κ B/p65-mediated caspase pathway. B CANCER. 2012 Apr 1;99(4):E49-54.
12. Meng X, Li F, Chen S, Tang C, Zhang W, Wang Z, Zhao S . Cloning and expression of neuron-specific enolase in the corpus luteum of dairy goats. Gene. 2012 Jul 25;503(2):222-8.
13. Liu Q, Wu C, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladder dysfunction associated with spinal cord injury. Int J Mol Med. 2018 May;41(5):2609-2618
14. Liu Q, Wu C, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladder dysfunction associated with spinal cord injury. Int J Mol Med. 2018 May;41(5):2609-2618.
15. Minghao Luo, Suxin Luo, Yuzhou Xue, Qing Chang, Hui Yang, Wenyu Dong, Ting Zhang, Shuyuan Cao . Aerobic exercise inhibits renal EMT by promoting irisin expression in SHR iScience. 2023 Jan 14;26(2):105990.
16. Xiumin Li, Qiaoling Su, Yutian Pan . Overcharged lipid metabolism in mechanisms of antitumor by Tremella fuciformis-derived polysaccharide Int J Oncol. 2023 Jan;62(1):11.
17. Yuming Zhang, Hai'ou Lu, Yuexin Yu . ZAP70 interaction with 13 mRNAs as a potential immunotherapeutic target for endometrial cancer Oncol Lett. 2023 Apr 6;25(5):213.
18. Huiling Zhang, Xuedi Zhang, Chunxiu Ling, Change Liu, Shaodong Hua, Ziyang Xiong, Huimin Zhao, Qiushuang Feng, Zhichun Feng, Jing Tang, Xiaoyang Hong . EGFR-TNFR1 pathway in endothelial cell facilitates acute lung injury by NF- κ B/MAPK-mediated inflammation and RIP3-dependent necroptosis Int Immunopharmacol. 2023 Apr;117:109902.
19. Lijuan Chen, Jing Li, Xinglian Liu, Zhiwei Zhao, Yan Jin, Yikun Fu, Aiqin Zhou, Chengqun Wang, Yan Zhou . Vitamin B6 Deficiency Induces Autism-Like Behaviors in Rats by Regulating mTOR-Mediated Autophagy in the Hippocampus Behav Neurol. 2023 May 9;2023:6991826.

Version 2025.03.12