BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)

产品编号	产品名称	包装
P2151-200μl	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	200μl
P2151-1ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	1ml
P2151-5ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	5ml

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Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer
Beads size	~200nm
Magnetization	Superparamagnetic
Coupled protein	Streptavidin
M.W. of protein	~55kDa (Streptavidin)
Protein concentration	≥ 0.6mg streptavidin per ml beads
Binding capacity	per mg beads: ≥20µg biotinylated antibody or dsDNA, ≥1000pmol free Biotin, ≥400pmol biotinylated oligonucleotides or peptides
Specificity	Biotinylated ligands
Elution method	Elution with acid or SDS-PAGE loading buffer for single-use applications
Application	purification of biotin-labeled proteins and nucleic acids, IP, Co-IP, DNA-protein pulldowns

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TBS	Tris Buffered Saline (e.g. ST661/ST665)	Prewash
Binding & Washing Buffer I (2X)	10mM Tris-HCl (pH 7.5), 1mM EDTA, 2M NaCl, 0.01%-0.1% Tween-20	for nucleic acid
Washing Buffer II (1X)	PBS (pH7.4), 0.05% Tween-20, with or without 0.01%-0.1% BSA	for antibody/protein
DNA/RNA Elution Buffer	95% formamide, 10mM EDTA, pH 8.2	for nucleic acid
Acid Elution Buffer	0.1M Glycine (pH2.0)	for antibody/protein
SDS-PAGE Loading Buffer	SDS-PAGE Loading Buffer (e.g. P0015)	for antibody/protein

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P2105-5ml	BeyoMag™ Protein G磁珠	5ml
P2108-1ml	BeyoMag™ Protein A+G磁珠	1ml
P2108-5ml	BeyoMag™ Protein A+G磁珠	5ml
P2115-0.5ml	BeyoMag™ Anti-Flag Magnetic Beads (Anti-Flag磁珠)	0.5ml
P2115-2ml	BeyoMag™ Anti-Flag Magnetic Beads (Anti-Flag磁珠)	2ml
P2118-0.5ml	BeyoMag™ Anti-Myc Magnetic Beads (Anti-Myc磁珠)	0.5ml
P2118-2ml	BeyoMag™ Anti-Myc Magnetic Beads (Anti-Myc磁珠)	2ml
P2121-0.5ml	BeyoMag™ Anti-HA Magnetic Beads (Anti-HA磁珠)	0.5ml
P2121-2ml	BeyoMag™ Anti-HA Magnetic Beads (Anti-HA磁珠)	2ml

Version 2021.03.01