碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Yun-Long Wang, Wan-Wen Zhao, Shao-Mei Bai, Li-Li Feng, Shu-Ying Bie, Li Gong, Fang Wang, Ming-Biao Wei, Wei-Xing Feng, Xiao-Lin Pang, Cao-Litao Qin, Xin-Ke Yin, Ying-Nai Wang, Weihua Zhou, Daniel R Wahl, Quentin Liu, Ming Chen, Mien-Chie Hung, Xiang-Bo Wan. MRNIP condensates promote DNA double-strand break sensing and end resection. *Nat Commun.* 2022 May 12;13(1):2638.
2. Niu K, Xiang L, Jin Y, Peng Y, Wu F, Tang W, Zhang X, Deng H, Xiang H, Li S, Wang J, Song Q, Feng Q. Identification of LARK as a novel and conserved G-quadruplex binding protein in invertebrates and vertebrates. *Nucleic Acids Res.* 2019 Aug 22;47(14):7306-7320.
3. Haitao Wang, Jianhua Zhang, Haoqiu Liu, Man Wang, Yan Dong, Yijun Zhou, Sek-Man Wong, Kai Xu, Qiufang Xu. A plant virus hijacks phosphatidylinositol-3,5-bisphosphate to escape autophagic degradation in its insect vector. *Autophagy.* 2023 Apr;19(4):1128-1143.
4. Jie Yao, Shuming Tang, Chenyan Shi, Yunzhi Lin, Lanlan Ge, Qinghua Chen, Baoru Ou, Dongyu Liu, Yuyang Miao, Qiuqie Xie, Xudong Tang, Jia Fei, Guangyi Yang, Jun Tian, Xiaobin Zeng. Isoginkgetin, a potential CDK6 inhibitor, suppresses SLC2A1/GLUT1 enhancer activity to induce AMPK-ULK1-mediated cytotoxic autophagy in hepatocellular carcinoma. *Autophagy.* 2023 Apr;19(4):1221-1238.
5. Liu Y, Long YH, Wang SQ, Zhang YY, Li YF, Mi JS, Yu CH, Li DY, Zhang JH, Zhang XJ. JMJD6 regulates histone H2A.X phosphorylation and promotes autophagy in triple-negative breast cancer cells via a novel tyrosine kinase activity. *Oncogene.* 2019 Feb;38(7):980-997.
6. Cheng-Wen He, Xue-Fei Cui, Shao-Jie Ma, Qin Xu, Yan-Peng Ran, Wei-Zhi Chen, Jun-Xi Mu, Hui Li, Jing Zhu, Qingqiu Gong, Zhiping Xie. Membrane recruitment of Atg8 by Hfl1 facilitates turnover of vacuolar membrane proteins in yeast cells approaching stationary phase. *BMC Biol.* 2021 Jun 4;19(1):117.
7. Tingting Zou, Yang Wang, Ling Dong, Tiantian Che, Huakan Zhao, Xiaohua Yan, Zhenghong Lin. Stabilization of SETD3 by deubiquitinase USP27 enhances cell proliferation and hepatocellular carcinoma progression. *Cell Mol Life Sci.* 2022 Jan 12;79(1):70.
8. Jin-Yu Liu, Ya-Jing Chen, Huan-Hui Feng, Zhan-Li Chen, Yun-Long Wang, Jin-E Yang, Shi-Mei Zhuang. LncRNA SNHG17 interacts with LRPPRC to stabilize c-Myc protein and promote G1/S transition and cell proliferation. *Cell Death Dis.* 2021 Oct 20;12(11):970.
9. Hanlin Ma, Gonghua Qi, Fang Han, Wei Lu, Jiali Peng, Rongrong Li, Shi Yan, Cunzhong Yuan, Beihua Kong. HMGB3 promotes PARP inhibitor resistance through interacting with PARP1 in ovarian cancer. *Cell Death Dis.* 2022 Mar 24;13(3):263.
10. Hanlin Ma, Gonghua Qi, Fang Han, Jiali Peng, Cunzhong Yuan, Beihua Kong. PBK drives PARP inhibitor resistance through the TRIM37/NFκB axis in ovarian cancer. *Exp Mol Med.* 2022 Jul;54(7):999-1010.
11. Jie Cheng, Peng Zhao, Lin Zhu, Fang Zhu, Zhiqiang Tian, Zhongjian Shen, Xiaoming Liu, Xiaoxia Liu. Corazonin signaling modulates the synthetic activity of male accessory gland in *Grapholita molesta*. *Int J Biol Macromol.* 2022 Sep 1:216:446-455.
12. Liu L, Zhou Y, Qu M, Qiu Y, Guo X, Zhang Y, Liu T, Yang J, Yang Q. Structural and biochemical insights into the catalytic mechanisms of two insect chitin deacetylases of the carbohydrate esterase 4 family. *J Biol Chem.* 2019 Apr 12;294(15):5774-5783.
13. Feng J, Li D, Liu L, Tang Y, Du R. Interaction of the small GTP-binding protein (Rab7) with β-actin in *Litopenaeus vannamei* and its role in white spot syndrome virus infection. *FISH SHELLFISH IMMUN.* 2019 May;88:1-8.
14. Feng J, Li D, Tang Y, Du R, Liu L. Molecular cloning of the Rab7 effector RILP (Rab-interacting lysosomal protein) in *Litopenaeus vannamei* and preliminary analysis of its role in white spot syndrome virus infection. *FISH SHELLFISH IMMUN.* 2019 Jul;90:126-133.
15. Deng S, Yang C, Nie K, Fan S, Zhu M, Zhu J, Chen Y, Yuan J, Zhang J, Xu H, Tian S, Chen J, Zhao M. Host cell protein PSMB10 interacts with viral NS3 protein and inhibits the growth of classical swine fever virus. *Virology.* 2019 Nov;537:74-83.
16. Ying-Ying Wu, Xiao-Fen Liu, Bei-Ling Fu, Qiu-Yun Zhang, Yang Tong, Jian Wang, Wen-Qiu Wang, Donald Grierson, Xue-Ren Yin. Methyl Jasmonate Enhances Ethylene Synthesis in Kiwifruit by Inducing NAC Genes That Activate ACS1. *J Agric Food Chem.* 2020 Mar 11;68(10):3267-3276.
17. Jianqiao Tu, Bijun Zhang, Guicun Fang, Wenjing Chang, Yanyan Zhao. Neddylation-mediated Nedd4-2 activation regulates ubiquitination modification of renal NBCe1. *Exp Cell Res.* 2020 May 15;390(2):111958.
18. Zijin Liu, Lamei Zheng, Li Pu, Xiaofeng Ma, Xing Wang, Yu Wu, Hainan Ming, Qing Wang, Genfa Zhang. ENO2 Affects the Seed Size and Weight by Adjusting Cytokinin Content and Forming ENO2-bZIP75 Complex in *Arabidopsis thaliana*. *Front Plant Sci.* 2020 Aug 26;11:574316.
19. Li-Zhen Zhang, Jin-E Yang, Yu-Wei Luo, Feng-Ting Liu, Yun-Fei Yuan, Shi-Mei Zhuang. A p53/Inc-Ip53 Negative Feedback Loop Regulates Tumor Growth and Chemoresistance. *Adv Sci (Weinh).* 2020 Sep 28;7(21):2001364.
20. Jindai Fan, Mengru Zhang, Chenchen Liu, Mengjiao Zhu, Zilin Zhang, Keke Wu, Zhaoyao Li, Wenhui Li, Shuangqi Fan, Chunmei Ju, Lin Yi, Hongxing Ding, Mingqiu Zhao, Jinding Chen. The Network of Interactions Between Classical Swine Fever Virus Nonstructural Protein p7 and Host Proteins. *Front Microbiol.* 2020 Nov 30;11:597893.
21. Hong Zhou, Fuyou Guo, Jinxiang Luo, Yongqiang Zhang, Jinlin Liu, Yanchun Zhang, Xinyu Zheng, Fenglin Wan, Wei Ding. Functional analysis of an upregulated calmodulin gene related to the acaricidal activity of curcumin against *Tetranychus cinnabarinus* (Boisduval). *Pest Manag Sci.* 2021 Feb;77(2):719-730.
22. Peifu Feng, Ling Li, Jing Dai, Lingli Zhou, Jing Liu, Jinfeng Zhao, Xiaodong Li, Neng Ling, Siyuan Qiu, Lin Zhang, Tiantian Xie, Yinglei Chen, Michael J Donovan, Tianhuan Peng, Jianhui Song, Mao Ye. The regulation of NONO by USP11 via deubiquitination is linked to the proliferation of melanoma cells. *J Cell Mol Med.* 2021 Feb;25(3):1507-1517.
23. Yu Huang, Qi Zheng, Zhiwen Wang, Jufen Tang, Yishan Lu, Qiwei Qin, Jia Cai, Jichang Jian. Fish natural killer enhancing factor-A (NKEF-A) enhance cytotoxicity of nonspecific cytotoxic cells against bacterial infection. *Mol Immunol.* 2021 May;133:184-193.
24. Dongfang Wang, Keai Sinn Tan, Xabier Arias-Moreno, Wen Tan, Guohua Cheng. Sirt3 increases CNPase enzymatic activity through deacetylation and facilitating substrate accessibility. *Biochem Biophys Res Commun.* 2021 Sep 24;571:181-187.
25. Yan Hu, Ya-Ping Liang, Hang-Yu Tian, Cui-Xian Xu, Dan Yu, Pan Zhang, Hui Ye, Ming Li. *Microplitis bicoloratus* bracovirus regulates cyclophilin A-apoptosis-inducing factor interaction to induce cell apoptosis in the insect immunosuppressive process. *Arch Insect Biochem Physiol.* 2022 May;110(1):e21877.
26. Xiongnan Chen, Xiaojun Chen, Yifan Liang, Sijia Xu, Zhijun Weng, Qi Gao, Zhao Huang, Guihong Zhang, Lang Gong. Interaction network of African swine fever virus structural protein p30 with host proteins. *Front Microbiol.* 2022 Aug 15:13:971888.
27. Luping Gong, Shenghao Liao, Wen Duan, Yongqiang Liu, Dongmei Zhu, Xiaosheng Zhou, Baoping Xue, Chengcai Chu, Yun-Kuan Liang. OsCPL3 is involved in brassinosteroid signaling by regulating OsGSK2 stability. *J Integr Plant Biol.* 2022 Aug;64(8):1560-1574.
28. Shan Sun, Hongyu Hou, Guoqiang Ma, Qilian Ma, Ningning Li, Li Zhang, Chunsheng Dong, Mian Cao, Kin Yip Tam, Zheng Ying, Hongfeng Wang. The interaction between E3 ubiquitin ligase Parkin and mitophagy receptor PHB2 links inner mitochondrial membrane ubiquitination to efficient mitophagy. *J Biol Chem.* 2022 Dec;298(12):102704.