碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)

产品编号	产品名称	包装
P2282-0.5ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	0.5ml
P2282-2ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	2ml
P2282-10ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	10ml

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Cat. No.	P2282	P2271
Product name	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)
Product content	50% settled gel in 50% glycerol with 10 mM phosphate-buffered saline and preservative (pH7.4)	50% settled gel in TBS with preservative (pH7.6)
Matrix	4% agarose	4% cross-linked agarose
Average bead size	~90µm	
Antibody	Mouse monoclonal antibody against Flag-tag	
Isotype	IgG2b	
M.W. of antibody	Approximately 150kDa	
Antibody concentration	Approximately 8mg Flag antibody per ml settled gel	
Binding capacity	Approximately 1mg Flag-tagged protein per ml settled gel	
Elution method	Acid, alkaline, neutral, peptide competitive or SDS-PAGE loading buffer elution. Note: If elute with SDS-PAGE loading buffer, the light (~25kDa) and heavy (~50kDa) chain of antibody will be denatured and release from the gel.	
Reagent	Chaotropic reagents will denature the target Flag-tagged protein. Do not exceed 0.3M GuHCl or 1.5M	

compatibility	Urea.	
Application	Suitable for IP、Co-IP and protein purification.	
Storage	-20°C	4°C

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Problem	Possible Causes	Solution
Large amount of tagged protein found in the flow through.	Binding time is not enough.	If using batch method, increase the binding time experimentally; If using column method, use a lower flow rate when loading samples.
	Column is overloaded.	Reduce the amount of the sample added to the gel or increase the amount of gel.
	Tag is not accessible to gel.	Expose the epitope tag by adding low amount of denaturant to the protein extract (dialysis may be needed before applying onto gel), or fuse the tag to the other terminus of the target protein.
	Gel has not been regenerated since last purification.	Perform gel regeneration procedure prior to binding.
	Reagent compatibility problem.	Dialyze the sample against TBS before purification procedure.
	The target protein has been degraded.	1. Prepare fresh lysates. Avoid using frozen lysates. 2. Perform purification at lower temperature, such as 4°C. 3. Use appropriate protease inhibitors in the lysate or increase their concentrations to prevent degradation of the fusion protein.
Very few or no tagged protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the tag by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	1. Reduce the number of washes. 2. Avoid adding high concentrations of NaCl to the mixture. 3. Use solutions that contain less or no detergent
	Incubation times are inadequate.	Increase the incubation times with the affinity gel (from several hours to overnight).
	Interfering substance is present in sample.	1. Lysates containing high concentrations of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided. 2. Excessive detergent concentrations may interfere with the antibody-antigen interaction. Detergent levels in buffers may be reduced by dilution.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity.

		2. Verify that the transfer was adequate by staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non-specific binding.	1. Prepare cell lysate again. 2. Add additional wash steps.
Background is too high.	Proteins bind nonspecifically to the monoclonal antibody, the gel beads, or the microcentrifuge tubes.	1. Pre-clear lysate with Mouse IgG-Agarose (P2055+A7028) to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before centrifugation.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins from the lysate during the initial centrifugation of the affinity gel complexes.

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产品编号	产品名称	包装
P2006	Protein A Agarose (Fast Flow, 进口分装)	2ml
P2009	Protein G Agarose (Fast Flow, 进口分装)	2ml
P2012	Protein A+G Agarose (Fast Flow, 进口分装)	2ml
P2051-2ml	Protein A Agarose (Fast Flow, for IP)	2ml
P2051-10ml	Protein A Agarose (Fast Flow, for IP)	10ml
P2051-50ml	Protein A Agarose (Fast Flow, for IP)	50ml
P2053-2ml	Protein G Agarose (Fast Flow, for IP)	2ml
P2053-10ml	Protein G Agarose (Fast Flow, for IP)	10ml
P2053-50ml	Protein G Agarose (Fast Flow, for IP)	50ml
P2055-2ml	Protein A+G Agarose (Fast Flow, for IP)	2ml
P2055-10ml	Protein A+G Agarose (Fast Flow, for IP)	10ml
P2055-50ml	Protein A+G Agarose (Fast Flow, for IP)	50ml
P2265	Mouse IgG Agarose (小鼠IgG琼脂糖凝胶)	1ml/5ml
P2267	Rabbit IgG Agarose (兔IgG琼脂糖凝胶)	1ml/5ml
P2271-0.5ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	0.5ml
P2271-2ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	2ml
P2271-10ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	10ml
P2282-0.5ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	0.5ml
P2282-2ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	2ml
P2282-10ml	Anti-Flag Affinity Gel (Anti-Flag亲和凝胶)	10ml
P2285-0.5ml	Anti-Myc Affinity Gel (Anti-Myc亲和凝胶)	0.5ml
P2285-2ml	Anti-Myc Affinity Gel (Anti-Myc亲和凝胶)	2ml
P2285-10ml	Anti-Myc Affinity Gel (Anti-Myc亲和凝胶)	10ml
P2287-0.5ml	Anti-HA Affinity Gel (Anti-HA亲和凝胶)	0.5ml
P2287-2ml	Anti-HA Affinity Gel (Anti-HA亲和凝胶)	2ml
P2287-10ml	Anti-HA Affinity Gel (Anti-HA亲和凝胶)	10ml
P2289-0.5ml	Anti-V5 Affinity Gel (Anti-V5亲和凝胶)	0.5ml
P2289-2ml	Anti-V5 Affinity Gel (Anti-V5亲和凝胶)	2ml
P2289-10ml	Anti-V5 Affinity Gel (Anti-V5亲和凝胶)	10ml

Version 2021.07.21