碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

Monomeric Streptavidin (单体链霉亲和素, powder)

产品编号	产品名称	包装
P5089-1mg	Monomeric Streptavidin (单体链霉亲和素, powder)	1mg
P5089-5mg	Monomeric Streptavidin (单体链霉亲和素, powder)	5mg
P5089-20mg	Monomeric Streptavidin (单体链霉亲和素, powder)	20mg
P5089-100mg	Monomeric Streptavidin (单体链霉亲和素, powder)	100mg

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About this protein	
Name	Monomeric Streptavidin (Recombinant monoStreptavidin; mSA); 单体链霉亲和素
Synonyms	N/A
Purity	>98% by SDS-PAGE
Biological Activity	>8U/mg
A282	2.8 (0.1% solution)
Formulation	lyophilized powder
Endotoxin	Less than 0.1EU/μg of monomeric streptavidin as determined by LAL method.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to precipitate the powder to the bottom. Reconstitute in ultrapure water or PBS to a concentration of 1mg/ml. Stock solutions should be apportioned into working aliquots and stored at ≤-20°C. Further dilutions should be made in appropriate buffered solutions.
Category	Streptavidin
Background	Monomeric Streptavidin is a mutant of Streptavidin, which is engineered by introducing mutations at the subunit interface and redesigning the binding pocket. The molecular mass of monomeric streptavidin is ~15kDa. Monomeric Streptavidin has only one binding site for biotin, but its engineering dissociation constant (KD) is like that of streptavidin [1-4].
Amino Acid Sequence	MHHHHHHGSA EAGITGTWY NQRGSTFTVT AGADGNLTGQ YENRAQGWGC QNSPYTLTGR YNGTKLEWRV EWNNSTENCH SRTEWRGQY QGGAEARINT QWNLTYEKGS GPATTQGQDT FTKVKPSAAS GSDYKDDDDK

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产品编号	产品名称	包装
P5089-1mg	Monomeric Streptavidin (powder)	1mg

P5089-5mg	Monomeric Streptavidin (powder)	5mg
P5089-20mg	Monomeric Streptavidin (powder)	20mg
P5089-100mg	Monomeric Streptavidin (powder)	100mg
—	说明书	1份

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