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1. Zhou F, Zhang L, Gong K, Lu G, Sheng B, Wang A, Zhao N, Zhang X, Gong Y. LEF-1 activates the transcription of E2F1. Biochem Biophys Res Commun. 2008 Jan 4;365(1):149-53.
2. Chen H, Sun JG, Cao XW, Ma XG, Xu JP, Luo FK, Chen ZT. Preliminary validation of ERBB2 expression regulated by miR-548d-3p and miR-559. Biochem Biophys Res Commun. 2009 Aug 7; 385(4):596-600.
3. Tian YP, Jiang W, Gao N, Zhang JL, Chen W, Fan DY, Zhou DS, An J. Inhibitory effects of glutathione on dengue virus production. Biochemical and Biophysical Research Communications. 2010 Jul 2; 397(3):420-4.
4. Zhang J, Chen Y, Xin XL, Li QN, Li M, Lin LP, Geng MY, Ding J. Oligomannurate sulfate blocks tumor growth by inhibiting NF-kappaB activation. Acta Pharmacol Sin. 2010 Mar;31(3):375-81.
5. Zhang J, Xin X, Chen Q, Xie Z, Gui M, Chen Y, Lin L, Feng J, Li Q, Ding J, Geng M. Oligomannurate sulfate sensitizes cancer cells to doxorubicin by inhibiting atypical activation of NF-kB via targeting of Mre11. Int J Cancer. 2012 Jan 15;130(2):467-77.
6. Hou S, Shu Q, Jiang Z, Chen Y, Li F, Chen F, Kijlstra A, Yang P. Replication study confirms the association between UBAC2 and Behçet's disease in two independent Chinese sets of patients and controls. Arthritis Res Ther. 2012 Mar 29;14(2):R70.
7. Zhang J, Chen J, Yang J, Xu CW, Pu P, Ding JW, Jiang H. Resveratrol attenuates oxidative stress induced by balloon injury in the rat carotid artery through actions on theERK1/2 and NF-kappa B pathway. Cell Physiol Biochem. 2013;31(2-3):230-41.
8. Liu D, Duan X, Dong D, Bai C, Li X, Sun G, Li B. Activation of the nrf2 pathway by inorganic arsenic in human hepatocytes and the role of

- transcriptional repressorbach1. *Oxid Med Cell Longev.* 2013; 2013:984546.
9. Ding Q, Zhong H, Qi Y, Cheng Y, Li W, Yan S, Wang X. Anti-arthritis effects of crocin in interleukin-1 β -treated articular chondrocytes and cartilage in a rabbitosteoarthritic model. *Inflamm Res.* 2013 Jan; 62(1):17-25.
 10. Yang CG, Liu SS, Sun B, Wang XL, Wang N, Chen SL. Iron-metabolic function and potential antibacterial role of Hepcidin and its correlated genes (Ferroportin 1 and Transferrin Receptor) in turbot (*Scophthalmus maximus*). *Fish Shellfish Immunol.* 2013 Mar;34(3):744-55.
 11. Song S, Zhou F, Chen WR, Xing D. PDT-induced HSP70 externalization up-regulates NO production via TLR2 signal pathway in macrophages. *FEBS Lett.* 2013 Jan 16;587(2):128-35.
 12. Deng X, Rui W, Zhang F, Ding W. PM2.5 induces Nrf2-mediated defense mechanisms against oxidative stress by activating PI3K/AKT signaling pathway in human lung alveolar epithelial A549 cells. *Cell Biol Toxicol.* 2013 Jun;29(3):143-57.
 13. Li C, Pan Z, Xu T, Zhang C, Wu Q, Niu Y. Puerarin induces the upregulation of glutathione levels and nuclear translocation of Nrf2 through PI3K/Akt/GSK-3 β signaling events in PC12 cells exposed to lead. *Neurotoxicol Teratol.* 2014 Sep 4;46C:1-9.
 14. Gao X, He Y, Gao LM, Feng J, Xie Y, Liu X, Liu L. Ser9-phosphorylated GSK3 β induced by 14-3-3 ζ actively antagonizes cell apoptosis in a NF- κ B dependent manner. *Biochem Cell Biol.* 2014 Oct;92(5):349-56.
 15. Li H, Liu F, Guo H, Zhu Z, Jiao Y. Role of interferon-inducible protein 202 (p202) in the regulation of adipogenesis in mouse adipose-derived stem cells. *Mol Cell Endocrinol.* 2014 Feb 15;382(2):814-24.
 16. Ruan GX, Chen YZ, Yao XL, Du A, Tang GP, Shen YQ, Tabata Y, Gao JQ. Macrophage mannose receptor-specific gene delivery vehicle for macrophage engineering. *Acta Biomater.* 2014 May;10(5):1847-55.
 17. Zeng W, Chang H, Ma M, Li Y. CCL20/CCR6 promotes the invasion and migration of thyroid cancer cells via NF-kappa B signaling-induced MMP-3 production. *Exp Mol Pathol.* 2014 Aug;97(1):184-90.
 18. Huang KF, Huang XP, Xiao GQ, Yang HY, Lin JS, Diao Y. Kallistatin, a novel anti-angiogenesis agent, inhibits angiogenesis via inhibition of the NF- κ B signaling pathway. *Biomed Pharmacother.* 2014 May; 68(4):455-61.
 19. Qin W, Pan Y, Zheng X, Li D, Bu J, Xu C, Tang J, Cui R, Lin P, Yu X. MicroRNA-124 regulates TGF- α -induced epithelial-mesenchymal transition in human prostate cancer cells. *Int J Oncol.* 2014 Sep; 45(3):1225-31.
 20. Han L, Tang C, Yin C. Oral delivery of shRNA and siRNA via multifunctional polymeric nanoparticles for synergistic cancer therapy. *Biomaterials.* 2014 May;35(15):4589-600.
 21. Ruan GX, Zhang TY, Li LM, Zhang XG, Shen YQ, Tabata Y, Gao JQ. Hepatic-targeted gene delivery using cationic mannan vehicle. *Mol Pharm.* 2014 Oct 6;11(10):3322-9.
 22. Li K, Feng L, Shen J, Zhang Q, Liu Z, Lee ST, Liu J. Patterned substrates of nano-graphene oxide mediating highly localized and efficient gene delivery. *ACS Appl Mater Interfaces.* 2014 Apr 23;6(8):5900-7.
 23. Wu WY, Yan H, Wang XB, Gui YZ, Gao F, Tang XL, Qin YL, Su M, Chen T, Wang YP. Sodium tanshinone IIA silicate inhibits high glucose-induced vascular smooth muscle cell proliferation and migration through activation of AMP-activated protein kinase. *PLoS One.* 2014 Apr 16;9(4):e94957.
 24. Gao F, Wang W. MicroRNA-96 promotes the proliferation of colorectal cancer cells and targets tumor protein p53 inducible nuclear protein 1, forkhead box protein O1 (FOXO1) and FOXO3a. *Mol Med Rep.* 2015 Feb;11(2):1200-6.
 25. Cao G, Gong S, Zhang F, Fu W. Xiao Yao San against Corticosterone-Induced Stress Injury via Upregulating Glucocorticoid Receptor Reaction Element Transcriptional Activity. *Evid Based Complement Alternat Med.* 2016;2016:5850739.
 26. Ding L, Yang M, Zhao T, Lv G. Roles of p300 and cyclic adenosine monophosphate response element binding protein in high glucose-induced hypoxia-inducible factor 1 α inactivation under hypoxic conditions. *J Diabetes Investig.* 2016 Nov 3. doi: 10.1111/jdi.12592. [Epub ahead of print]
 27. Chen Z, Yang L, Cui Y, Zhou Y, Yin X, Guo J, Zhang G, Wang T, He QY. Cytoskeleton-centric protein transportation by exosomes transforms tumor-favorable macrophages. *Oncotarget.* 2016 Oct 11; 7(41):67387-67402.
 28. Liu J, Li Y, Luo M, Yuan Z, Liu J. MicroRNA-214 inhibits the osteogenic differentiation of human osteoblasts through the direct regulation of baculoviral IAP repeat-containing 7. *Exp Cell Res.* 2017 Jan 18. pii: S0014-4827(17)30016-2.

Version 2022.09.14