碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-168-3301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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- Kong Y, Wei ZF, Fu YJ, Gu CB, Zhao CJ, Yao XH, Efferth Thomas. Negative-pressure cavitation extraction of cajaninstilbene acid and pinostrobin from pigeon pea [*Cajanus cajan* (L.) Millsp.] leaves and evaluation of antioxidant activity. *Food Chemistry*. 2011 Oct;128(3):596-605.
- Li W, Wu Y, Ren C, Lu Y, Gao Y, Zheng X, Zhang C. The activity of recombinant human neuroglobin as an antioxidant and free radical scavenger. *Proteins*. 2011 Jan;79(1):115-25.
- Wang L, Hao J, Hu J, Pu J, Lv Z, Zhao L, Wang Q, Yu Q, Wang Y, Li G. Protective Effects of Ginsenosides against Bisphenol A-Induced Cytotoxicity in

- 15P-1 Sertoli Cells via Extracellular Signal-Regulated Kinase 1/2 Signalling and Antioxidant Mechanisms. *Basic Clin Pharmacol Toxicol*. 2012 Jul; 111(1):42-9.
4. Zhou Y, Yuan J, Wu J, Han X. The toxic effects of microcystin-LR on rat spermatogonia in vitro. *Toxicol Lett*. 2012 Jul 7;212(1):48-56.
  5. Luo JG, Li L, Kong LY. Preparative separation of phenylpropanoid glycerides from the bulbs of *Lilium lancifolium* by high-speed counter-current chromatography and evaluation of their antioxidant activities. *Food Chemistry*. 2012 Apr;131(1):1056-62.
  6. Liu F, Li XL, Lin T, He DW, Wei GH, Liu JH, Li LS. The cyclophosphamide metabolite, acrolein, induces cytoskeletal changes and oxidative stress in Sertoli cells. *Mol Biol Rep*. 2012 Jan;39(1):493-500.
  7. Yang T, Luo F, Shen Y, An J, Li X, Liu X, Ying B, Liao Z, Dong J, Guo L, Wang T, Xu D, Chen L, Wen F. Quercetin attenuates airway inflammation and mucus production induced by cigarette smoke in rats. *Int Immunopharmacol*. 2012 May;13(1):73-81.
  8. Zhang JX, Yang JR, Chen GX, Tang LJ, Li WX, Yang H, Kong X. Sesamin ameliorates arterial dysfunction in spontaneously hypertensive rats via downregulation of NADPH oxidase subunits and upregulation of eNOS expression. *Acta Pharmacol Sin*. 2013 Jul 5;34(7):912-20.
  9. Yang L, Qu M, Wang Y, Duan H, Chen P, Wang Y, Shi W, Danielson P, Zhou Q. Trichostatin A inhibits transforming growth factor- $\beta$ -induced reactive oxygen species accumulation and myofibroblast differentiation via enhanced NF-E2-related factor 2-antioxidant response element signaling. *Mol Pharmacol*. 2013 Mar;83(3):671-80.
  10. Zou Y, Wang R, Guo H, Dong M. Phytoestrogen  $\beta$ -Ecdysterone Protects PC12 Cells Against MPP+-Induced Neurotoxicity In Vitro: Involvement of PI3K-Nrf2-Regulated Pathway. *Toxicol Sci*. 2015 Sep;147(1):28-38.
  11. Zhao P, Duan L, Guo L, Dou LL, Dong X, Zhou P, Li P, Liu EH. Chemical and biological comparison of the fruit extracts of *Citrus wilsonii* Tanaka and *Citrus medica* L. *Food Chem*. 2015 Apr 15;173:54-60.
  12. Zhang DY, Wan Y, Xu JY, Wu GH, Li L, Yao XH. Ultrasound extraction of polysaccharides from mulberry leaves and their effect on enhancing antioxidant activity. *Carbohydr Polym*. 2016 Feb 10;137:473-9.
  13. Yao C, Li G, Qian Y, Cai M, Yin H, Xiao L, Tang W, Guo F, Shi B. Protection of Pentoxifylline against Testis Injury Induced by Intermittent Hypobaric Hypoxia. *Oxid Med Cell Longev*. 2016;2016:3406802.
  14. Xia Q, Wang L, Xu C, Mei J, Li Y. Effects of germination and high hydrostatic pressure processing on mineral elements, amino acids and antioxidants in vitro bioaccessibility, as well as starch digestibility in brown rice (*Oryza sativa* L.). *Food Chem*. 2017 Jan 1;214:533-42.

Version 2017.03.09