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
1. Ji H, Wang L, Bi H, Sun L, Cai B, Wang Y, Zhao J, Du Z. Mechanisms of lumbrokinase in protection of cerebral ischemia. *Eur J Pharmacol*. 2008 Aug 20;590(1-3):281-9.
2. Bi J, Wang XB, Chen L, Hao S, An LJ, Jiang B, Guo L. Catalpol protects mesencephalic neurons against MPTP induced neurotoxicity via attenuation of mitochondrial dysfunction and MAO-B activity. *Toxicol In Vitro*. 2008 Dec;22(8):1883-9.
3. Guo RW, Wang H, Gao P, Li MQ, Zeng CY, Yu Y, Chen JF, Song MB, Shi YK, Huang L. An essential role for stromal interaction molecule 1 in neointima formation following arterial injury. *Cardiovasc Res*. 2009 Mar 1;81(4):660-8.
4. Yang LX, Guo RW, Liu B, Wang XM, Qi F, Guo CM, Shi YK, Wang H. Role of TRPC1 and NF-kappaB in mediating angiotensin II-induced Ca²⁺ entry and endothelial hyperpermeability. *Peptides*. 2009 Jul;30(7):1368-73.
5. Sun ZL, Niu RY, Su K, Wang B, Wang JM, Zhang JH, Wang JD. Effects of sodium fluoride on hyperactivation and Ca²⁺ signaling pathway in sperm from mice: an in vivo study. *ARCHIVES OF TOXICOLOGY*. 2010;84(5):353-61.
6. Shen W, Song D, Wu J, Zhang W. Protective effect of a polysaccharide isolated from a cultivated *Cordyceps* mycelia on hydrogenperoxide-induced oxidative damage in PC12 cells. *Phytother Res*. 2011 May;25(5):675-80.
7. Li C, Yang Z, Li Z, Ma Y, Zhang L, Zheng C, Qiu W, Wu X, Wang X, Li H, Tang J, Qian M, Li D, Wang P, Luo J, Liu M.. Maslinic acid suppresses osteoclastogenesis and prevents ovariectomy-induced bone loss by regulating RANKL-mediated NF-κB and MAPK signaling pathways. *J Bone Miner Res*. 2011 Mar;26(3):644-56.
8. Sun C, Wang L, Yan J, Liu S. Calcium ameliorates obesity induced by high-fat diet and its potential correlation with p38 MAPK pathway. *Mol Biol Rep*. 2012 Feb;39(2):1755-63.
9. Liu H, Han F, Shi Y. Effect of Calreticulin on Ca(2+)/CaM KinaseIIα and Endoplasmic Reticulum Stress in Hippocampal in a Rat Model of Post-traumatic Stress Disorder. *Neurochem Res*. 2013 Jul;38(7):1407-14.
10. Zhou J, Yang C, Wang J, Sun P, Fan P, Tian K, Liu S, Xia C. Toxic effects of environment-friendly antifoulant nonivamide on *Phaeodactylum tricornutum*. *Environ Toxicol Chem*. 2013 Apr;32(4):802-9.
11. Wang H, Yin G, Yu CH, Wang Y, Sun Z. Inhibitory effect of sanguinarine on PKC-CPI-17 pathway mediating by muscarinic receptors in dispersed intestinal smooth muscle cells. *Res Vet Sci*. 2013 Dec;95(3):1125-33.
12. Zheng B, Wu L, Ma L, Liu S, Li L, Xie W, Li X. Telekin induces apoptosis associated with the mitochondria-mediated pathway in human hepatocellular carcinoma cells. *Biol Pharm Bull*. 2013;36(7):1118-25.
13. Zhong L, Lv L, Yang J, Liao X, Yu J, Wang R, Zhou P. Inhibitory effect of hydrogen sulfide on platelet aggregation and the underlying mechanisms. *J Cardiovasc Pharmacol*. 2014 Nov;64(5):481-7.
14. Ren H, Teng Y, Tan B, Zhang X, Jiang W, Liu M, Jiang W, Du B, Qian M. Toll-like receptor-triggered calcium mobilization protects mice against bacterial infection through extracellular ATP release. *Infect Immun*. 2014 Dec;82(12):5076-85.
15. Zhou R, Ding XL, Liu LM. Ryanodine receptor 2 contributes to hemorrhagic shock-induced bi-phasic vascular reactivity in rats. *Acta Pharmacol Sin*. 2014 Nov;35(11):1375-84.
16. Liu Z, Cai H, Zhu H, Toque H, Zhao N, Qiu C, Guan G, Dang Y, Wang J. Protein kinase RNA-like endoplasmic reticulum kinase (PERK)/calcineurin signaling is a novel pathway regulating intracellular calcium accumulation which might be involved in ventricular arrhythmias in diabetic cardiomyopathy. *Cell Signal*. 2014 Dec;26(12):2591-600.
17. Liu J, Zhang Y, Li Q, Zhuang Q, Zhu X, Pan L, Cheng M. An improved method for guinea pig airway smooth muscle cell culture and the effect of SPFF on intracellular calcium. *Mol Med Rep*. 2014 Sep;10(3):1309-14.
18. Zhang Y, Yin H, Zhao X, Wang W, Du Y, He A, Sun K. The promoting effects of alginate oligosaccharides on root development in *Oryza sativa* L. mediated by auxin signaling. *Carbohydr Polym*. 2014 Nov 26;113:446-54.
19. Liu Y, Shen S, Li Z, Jiang Y, Si J, Chang Q, Liu X, Pan R. Cajanin stilbene acid protects corticosterone-induced injury in PC12 cells by inhibiting oxidative and endoplasmic reticulum stress-mediated apoptosis. *Neurochem Int*. 2014 Dec;78:43-52.
20. Xu X, Chen D, Ye B, Zhong F, Chen G. Curcumin induces the apoptosis of non-small cell lung cancer cells through a calcium signaling pathway. *Int J Mol Med*. 2015 Jun;35(6):1610-6.
21. Song X, Zhao C, Dai C, Ren Y, An N, Wen H, Pan LI, Cheng M, Zhang Y. Suppression of the increasing level of acetylcholine-stimulated intracellular Ca²⁺ in guinea pig airway smooth muscle cells by mabuterol. *Biomed Rep*. 2015 Nov;3(6):778-86.
22. Li H, Feng L, Jiang W, Liu Y, Jiang J, Zhang Y, Wu P, Zhou X. Ca(2+) and caspases are involved in hydroxyl radical-induced apoptosis in erythrocytes of Jian carp (*Cyprinus carpio* var. Jian). *Fish Physiol Biochem*. 2015 Oct;41(5):1305-19.
23. Wen L, Han F, Shi Y. Changes in the glucocorticoid receptor and Ca²⁺/calreticulin-dependent signalling pathway in the medial prefrontal cortex of rats with post-traumatic stress disorder. *J Mol Neurosci*. 2015 May;56(1):24-34.
24. Cao Z, Liu D, Zhang Q, Sun X, Li Y. Aluminum Chloride Induces Osteoblasts Apoptosis via Disrupting Calcium Homeostasis and Activating Ca²⁺/CaMKII Signal Pathway. *Biol Trace Elem Res*. 2016 Feb;169(2):547-53.
25. Wei NN, Lv HN, Wu Y, Yang SL, Sun XY, Lai R, Jiang Y, Wang K. Selective Activation of Nociceptor TRPV1 Channel and Reversal of Inflammatory Pain in Mice by a Novel Coumarin Derivative Muralatin L from *Murraya alata*. *J Biol Chem*. 2016 Jan 8;291(2):640-51.
26. Wang B, Zhao XH. Apigenin induces both intrinsic and extrinsic pathways of apoptosis in human colon carcinoma HCT-116 cells. *Oncol Rep*. 2017 Feb;37(2):1132-1140.

Version 2017.03.09