AAV Purification Kit Mini for All Serotypes

产品编号	产品名称	包装
C2906S	AAV Purification Kit Mini for All Serotypes	5次
C2906M	AAV Purification Kit Mini for All Serotypes	20次

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Table 1. Summary of Recommended Biosafety Levels for Infectious Agents

BSL	Agents	Practices	Primary Barriers and Safety Equipment	Facilities (Secondary Barriers)
1	Not known to consistently cause diseases in healthy adults	Standard microbiological practices	<ul style="list-style-type: none"> ■ No primary barriers required. ■ PPE: laboratory coats and gloves; eye, face protection, as needed 	Laboratory bench and sink required
2	<ul style="list-style-type: none"> ■ Agents associated with human disease ■ Routes of transmission include percutaneous injury, ingestion, mucous membrane exposure 	BSL-1 practice plus: <ul style="list-style-type: none"> ■ Limited access ■ Biohazard warning signs ■ ‘Sharps’ precautions ■ Biosafety manual defining any needed waste decontamination or medical surveillance policies 	Primary barriers: <ul style="list-style-type: none"> ■ BSCs or other physical containment devices used for all manipulations of agents that cause splashes or aerosols of infectious materials ■ PPE: Laboratory coats, gloves, face and eye protection, as needed 	BSL-1 plus: <ul style="list-style-type: none"> ■ Autoclave available
3	Indigenous or exotic agents that may cause serious or	BSL-2 practice plus: <ul style="list-style-type: none"> ■ Controlled access 	Primary barriers: <ul style="list-style-type: none"> ■ BSCs or other 	BSL-2 plus: <ul style="list-style-type: none"> ■ Physical separation

	potentially lethal disease through the inhalation route of exposure	<ul style="list-style-type: none"> ■ Decontamination of all waste ■ Decontamination of laboratory clothing before laundering 	physical containment devices used for all open manipulations of agents <ul style="list-style-type: none"> ■ PPE: Protective laboratory clothing, gloves, face, eye and respiratory protection, as needed 	from access corridors <ul style="list-style-type: none"> ■ Self-closing, double-door access ■ Exhausted air not recirculated ■ Negative airflow into laboratory ■ Entry through airlock or anteroom ■ Hand washing sink near laboratory exit
4	<ul style="list-style-type: none"> ■ Dangerous/exotic agents which post high individual risk of aerosol-transmitted laboratory infections that are frequently fatal, for which there are no vaccines or treatments ■ Agents with a close or identical antigenic relationship to an agent requiring BSL-4 until data are available to redesignate the level ■ Related agents with unknown risk of transmission 	BSL-3 practices plus: <ul style="list-style-type: none"> ■ Clothing change before entering ■ Shower on exit ■ All material decontaminated on exit from facility 	Primary barriers: <ul style="list-style-type: none"> ■ All procedures conducted in Class III BSCs or Class I or II BSCs in combination with full-body, air-supplied, positive pressure suit 	BSL-3 plus: <ul style="list-style-type: none"> ■ Separate building or isolated zone ■ Dedicated supply and exhaust, vacuum, and decontamination systems ■ Other requirements outlined in the text

BSL, biosafety level; PPE, personal protective equipment.

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产品编号	产品名称	包装
C2901	通用型病毒浓缩试剂盒	150ml/600ml/3L
C2906	AAV Purification Kit Mini for All Serotypes	5次/20次
C2908	AAV Purification Kit Midi for All Serotypes	5次/20次
C3013	AAV-mCherry-GFP-LC3B (AAV9, 10 ¹³ vp/ml)	250μl
C4107	AAV-CAG-GFP (AAV9, 10 ¹³ vp/ml)	250μl
C4108	AAV-CAG-Cre-GFP (AAV9, 10 ¹³ vp/ml)	250μl
FCD150	BeyoGold™ 150mm细胞培养皿	5个/袋, 12袋/箱
FFLK176	BeyoGold™ 175cm ² 透气盖细胞培养瓶	5个/包, 40个/箱
FSCP023	BeyoGold™ 23cm细胞刮(独立纸塑包装, 无菌)	100个/盒
FSCP029	BeyoGold™ 29cm细胞刮(独立纸塑包装, 无菌)	100个/盒
C3013	AAV-mCherry-GFP-LC3B (AAV9, 10 ¹³ vp/ml)	250μl
FUF058	超滤管(0.5ml, 100kDa MWCO, PES, Sartorius分装)	5个/25个/100个
FUF158	超滤管(15ml, 100kDa MWCO, PES, Sartorius分装)	2个/12个/48个
ST876	BeyoPure™ Ultrapure Water (DNase/RNase-Free, Sterile)	100ml/500ml

Version 2025.02.07