碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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1. Zhao Y, Chen X, Cai L, Yang Y, Sui G, Wu J. Angiotensin II suppresses adriamycin-induced apoptosis through activation of phosphatidylinositol 3-kinase/Akt signaling in human breast cancer cells. *Acta Biochim Biophys Sin (Shanghai)*. 2008 Apr;40(4):304-10.
2. Zhan X, Hu H, Ke C, Hu S, Wang D, Chen F. Isolation and characterization of eleven microsatellite loci in small abalone, *Haliotis diversicolor Reeve*. *Conserv Genet*. (2009) 10:1185-1187.
3. Shi G, Zhang Z, Feng D, Xu Y, Lu Y, Wang J, Jiang J, Zhang Z, Li X, Ning G. Selection of reference genes for quantitative real-time reverse transcription-polymerase chain reaction in concanavalin A-induced hepatitis model. *Anal Biochem*. 2010 Jun 1;401(1):81-90.
4. Ma L, Zhang XX, Cheng S, Zhang Z, Shi P, Liu B, Wu B, Zhang Y. Occurrence, abundance and elimination of class I integrons in one municipal sewage treatment plant. *Ecotoxicology*. 2011 Jul;20(5):968-73.
5. Shi G, Zhang Z, Zhang R, Zhang X, Lu Y, Yang J, Zhang D, Zhang Z, Li X, Ning G. Protective effect of andrographolide against concanavalin A-induced liver injury. *Naunyn Schmiedebergs Arch Pharmacol*. 2012 Jan;385(1):69-79.
6. Wan Z, Lu Y, Liao Q, Wu Y, Chen X. Fangchinoline Inhibits Human Immunodeficiency Virus Type 1 Replication by Interfering with gp160 Proteolytic Processing. *PLoS One*. 2012;7(6):e39225.
7. Shi H, Yan X, Ruan L, Xu X. A novel JNK from *Litopenaeus vannamei* involved in white spot syndrome virus infection. *Dev Comp Immunol*. 2012 Jul;37(3-4):421-8.
8. Zhou BR, Xu Y, Wu D, Permatasari F, Gao YY, Luo D. Ginsenoside Rg1 protects human fibroblasts against psoralen- and UVA-induced premature senescence through a telomeric mechanism. *Arch Dermatol Res*. 2012 Apr;304(3):223-8.
9. Liao R, Shen K, Li AM, Shi P, Li Y, Shi Q, Wang Z. High-nitrate wastewater treatment in an expanded granular sludge bed reactor and microbial diversity using 454pyrosequencing analysis. *Bioresour Technol*. 2013 Apr;134:190-7.
10. Hao G, Shi YH, Tang YL, Le GW. The intracellular mechanism of action on *Escherichia coli* of BF2-A/C, two analogues of the antimicrobial peptide Buforin 2. *J Microbiol*. 2013 Apr;51(2):200-6.
11. Xin T, Jia H, Ding J, Li P, Yang H, Hou S, Yuan W, Guo X, Wang H, Liang Q, Li M, Wang B, Zhu H. Assessment of a protein cocktail-based skin test for bovine tuberculosis in a double-blind field test in cattle. *Clin Vaccine Immunol*. 2013 Apr;20(4):482-90.
12. Dong B, Dai G, Xu L, Zhang Y, Ling L, Sun L, Lv J. Tumor cell lysate induces the immunosuppression and apoptosis of mouse immunocytes. *Mol Med Rep*. 2014 Dec;10(6):2827-34.
13. Shi Y, Wang D, Lu L, Yin Y, Wang M, Li C, Diao J, Wang Y, Wei L. Ligustilide prevents the apoptosis effects of tumour necrosis factor-alpha during C2C12 cell differentiation. *Int Immunopharmacol*. 2014 Apr;19(2):358-64.
14. He H, Xiao H, Kuang H, Xie Z, Chen X, Jing X, Huang Y. Synthesis of mesoporous silica nanoparticle-oxaliplatin conjugates for improved anticancer drug delivery. *Colloids Surf B Biointerfaces*. 2014 May 1; 117:75-81.
15. Chen LL, Han WF, Geng Y, Su JS. A genome-wide study of DNA methylation modified by epigallocatechin-3-gallate in the CAL-27 cell line. *Mol Med Rep*. 2015 Oct;12(4):5886-90.
16. Zhang Y, Bo Q, Wu W, Xu C, Yu G, Ma S, Yang Q, Cao Y, Han Q, Ru Y, Liu X, Hua Wei R, Wang FE, Zhang X, Li X.  $\alpha$ -Melanocyte-stimulating hormone prevents glutamate excitotoxicity in developing chicken retina via MC4R-mediated down-regulation of microRNA-194. *Sci Rep*. 2015 Oct 28;5:15812.
17. Ma J, Liu J, Lu C, Cai D. Pachymic acid induces apoptosis via activating ROS-dependent JNK and ER stress pathways in lung cancer cells. *Cancer Cell Int*. 2015 Aug 5;15:78.
18. Hu GL, Yan G, Xu H, Hua BZ. Molecular phylogeny of Panorpidae (Insecta: Mecoptera) based on mitochondrial and nuclear genes. *Mol Phylogenet Evol*. 2015 Apr;85:22-31.
19. Shi M, Zheng J, Tan Y, Tan G, Li J, Li Y, Li X, Zhou Z, Yang R. Ultrasensitive detection of single nucleotide polymorphism in human mitochondrial DNA utilizing ion-mediated cascade surface-enhanced Raman spectroscopy amplification. *Anal Chem*. 2015 Mar 3;87(5):2734-40.
20. Yang J, Yang C, Liu C, Zhang T, Yang Z. Paradoxical effects of VEGF on synaptic activity partially involved in notch1 signaling in the mouse hippocampus. *Hippocampus*. 2016 May;26(5):589-600.
21. Zhang X, Zhang P, Gao J, Huang Q. Autophagy dysregulation caused by ApoM deficiency plays an important role in liver lipid metabolic disorder. *BIOCHEM BIOPH RES CO*. 2018 Jan 22;495(4):2643-2648.
22. Zhang X, Zhang P, Gao J, Huang Q. Autophagy dysregulation caused by ApoM deficiency plays an important role in liver lipid metabolic disorder. *BIOCHEM BIOPH RES CO*. 2018 Jan 22;495(4):2643-2648.
23. Yang Z, Sun Y, Xian L, Xun Z, Yu J, Yang T, Zhao X, Cai C, Wang D, Ding P. Disulfide-bond-containing agmatine-cystaminebisacrylamide polymer demonstrates better transfection efficiency and lower cytotoxicity than polyethylenimine in NIH/3T3 cells. *J Cell Biochem*. 2018 Feb;119(2):1767-1779.
24. Han Y, Yin W, Li J, Zhao H, Zha Z, Ke W, Wang Y, He C, Ge Z. Intracellular glutathione-depleting polymeric micelles for cisplatin prodrug delivery to overcome cisplatin resistance of cancers. *J Control Release*. 2018 Mar 10;273:30-39.
25. Wang C, Wang F, Li Z, Cao Q, Huang L, Chen S. MeCP2-mediated epigenetic regulation in senescent endothelial progenitor cells. *Stem Cell Res Ther*. 2018 Apr 3;9(1):87.
26. Wang C, Wang F, Li Z, Cao Q, Huang L, Chen S. MeCP2-mediated epigenetic regulation in senescent endothelial progenitor cells. *Stem Cell Res Ther*. 2018 Apr 3.
27. Du X, Shen T, Wang H, Qin X, Xing D, Ye Q, Shi Z, Fang Z, Zhu Y, Yang Y, Peng Z, Zhao C, Lv B, Li X, Liu G, Li X. Adaptations of hepatic lipid metabolism and mitochondria in dairy cows with mild fatty liver. *J Dairy Sci*. 2018 Oct;101(10):9544-9558.
28. Wang W, Guo X, Li YM, Wang XY, Yang XJ, Wang YF, Wang TY. Enhanced transgene expression using cis-acting elements combined with the EF1 promoter in a mammalian expression system. *Eur J Pharm Sci*. 2018 Oct 15;123:539-545.
29. Guo Q, Zheng M, Xu Y, Wang N, Zhao W. MiR-384 induces apoptosis and autophagy of non-small cell lung cancer cells through the negative regulation of Collagen  $\alpha$ -1(X) chain gene. *BIOSCIENCE REP*. 2019 Feb 1;39(2). pii: BSR20181523.
30. Zhang J, Tao X, Sun M, Ying R, Su W, Wei W, Meng X. A Rat Model of Radiation Vasculitis for the Study of Mesenchymal Stem Cell-Based Therapy. *Biomed Res Int*. 2019 Mar 6;2019:3727635.
31. Ha S, Zhou H, Gautam M, Song Y, Wang C. Reduced ribosomal RNA expression and unchanged ribosomal DNA promoter methylation in oral squamous cell carcinoma. *Mol Genet Genomic Med*. 2019 Jul;7(7):e00783.
32. Fang Y, Zhu L, An N, Jiang G, Qian J, Zhao R, Yuan N, Zhang S, Wang J. Blood autophagy defect causes accelerated non-hematopoietic organ aging. *AGING-US*. 2019 Jul 21;11(14):4910-4922.
33. Zha S, Li Z, Chen S, Liu F, Wang F. MeCP2 inhibits cell functionality through FoxO3a and autophagy in endothelial progenitor cells. *AGING-US*. 2019 Sep 2;11(17):6714-6733.
34. Zheng Q, Zhang X, Yang H, Xie J, Xie Y, Chen J, Yu C, Zhong C. Internal Ribosome Entry Site Dramatically Reduces Transgene Expression in Hematopoietic Cells in a Position-Dependent Manner. *VIRUSES-BASEL*. 2019 Oct 8;11(10). pii: E920.
35. Zhu Y, Xing L, Zheng X, Yang CX, He YJ, Zhou TJ, Jin QR, Jiang HL. Amplification of tumor antigen presentation by NLGp1atin to improve chemioimmunotherapy. *INT J PHARMACOL*. 2020 Jan 5;573:118736.
36. Meng-Meng Ren, Sen Xu, Yu-Bo Wei, Juan-Juan Yang, Ya-Nan Yang, Shan-