



碧云天生物技术/Beyotime Biotechnology
订货热线：400-168-3301或800-8283301
订货e-mail：order@beyotime.com
技术咨询：info@beyotime.com
网址：http://www.beyotime.com

L-NAME (eNOS抑制剂)

产品编号	产品名称	包装
S0006	L-NAME (eNOS抑制剂)	200mg

产品简介：

- L-NAME, 即N^G-Nitro-L-arginine Methyl Ester, Hydrochloride, 也称L-NAME, HCl或N^ω-NO₂-L-Arg-OMe或N^G-NO₂-L-Arg-OMe, 是eNOS (endothelial nitric oxide synthase)的可逆抑制剂, IC₅₀=500nM。高浓度L-NAME(例如100μM或更高浓度)也可以显著抑制iNOS和nNOS。
- L-NAME为接近白色的粉末, 分子量269.7, 分子式为C₇H₁₅N₅O₄·HCl, 纯度大于98%。
- 溶解于水, 可以配制成50mg/ml的水溶液。

包装清单：

产品编号	产品名称	包装
S0006	L-NAME (eNOS抑制剂)	200mg
—	说明书	1份

保存条件：

-20°C避光保存, 一年有效。

注意事项：

- 如果配制成水溶液, 分装后-20°C保存, 4个月内有效。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明：

1. L-NAME储备液的配制: 称取少量L-NAME, 推荐使用Milli-Q级纯水溶解L-NAME, 使浓度为50mg/ml。
2. L-NAME的工作浓度通常为10-100μM。少数情况下浓度需达到500μM甚至5mM才可以观察到显著效果。其最佳工作浓度需根据具体的实验, 自行摸索。可以先分别尝试10、30和100μM这三个浓度。如果效果不佳, 可以尝试1μM等低浓度或300μM等更高浓度。

使用本产品的文献：

1. Li XL, Zou XM, Gao P, Li YL, Wang H, Chen XW. Role of nitric oxide in ischemia-reperfusion injury and acute rejection in rat intestinal transplantation. Transplant Proc. 2008 Dec;40(10):3342-5.
2. Wang GF, Wu SY, Rao JJ, Lü L, Xu W, Pang JX, Liu ZQ, Wu SG, Zhang JJ. Genipin inhibits endothelial exocytosis via nitric oxide in cultured human umbilical vein endothelial cells. Acta Pharmacol Sin. 2009 May;30(5):589-96.
3. Zhong W, Zou G, Gu J, Zhang J. L-arginine attenuates high glucose-accelerated senescence in human umbilical vein endothelial cells. Diabetes Res Clin Pract. 2010;89(1):38-45.
4. Zhang QH, Zu XY, Cao RX, Liu JH, Mo ZC, Zeng Y, Li YB, Xiong SL, Liu X, Liao DF, Yi GH. An involvement of SR-B1 mediated PI3K-Akt-eNOS signaling in HDL-induced cyclooxygenase 2 expression and prostacyclin production in endothelial cells. Biochem Biophys Res Commun. 2012 Mar 30;420(1):17-23.
5. Li ZL, Liu JC, Liu SB, Li XQ, Yi DH, Zhao MG Improvement of vascular function by acute and chronic treatment with the GPR30 agonist G1 in experimental diabetes mellitus. PLoS One. 2012;7(6):e38787.
6. Hong Q, Qi K, Feng Z, Huang Z, Cui S, Wang L, Fu B, Ding R, Yang J, Chen X, Wu D Hyperuricemia induces endothelial dysfunction via mitochondrial Na⁺/Ca²⁺ exchanger-mediated mitochondrial calcium overload. Cell Calcium. 2012 May;51(5):402-10.
7. Zhang Y, Liu H, Yin H, Wang W, Zhao X, Du Y. Nitric oxide mediates alginate oligosaccharides-induced root development in wheat (Triticum aestivum L.) Plant Physiology and Biochemistry.2013 Oct;71:49-56.
8. Guo Y, Guo M, Zhao W, Chen K, Zhang P. Burdock fructooligosaccharide induces stomatal closure in Pisum sativum Carbohydrate Polymers. 2013 Sep 12;97(2):731-735.
9. Chen W, Wang J, Chen J, Chen J, Chen Z. Relationship between changes in the cochlear blood flow and disorder of hearing function induced by blast injury in guinea pigs. Int J Clin Exp Pathol. 2013;6(3):375-84.
10. Ma S, Zhang Z, Yi F, Wang Y, Zhang X, Li X, Yuan Y, Cao F. Protective effects of low-frequency magnetic fields on cardiomyocytes from ischemia-reperfusion injury via ROS and NO/ONOO-. Oxid Med Cell Longev. 2013;2013:529173.
11. Pei H, Yu Q, Xue Q, Guo Y, Sun L, Hong Z, Han H, Gao E, Qu Y, Tao L. Notch1 cardioprotection in myocardial ischemia/reperfusion involves reduction of oxidative/nitrative stress. Basic Res Cardiol. 2013 Sep;108(5):373.
12. Yan T, Cui K, Huang X, Ding S, Zheng Y, Luo Q, Liu X, Zou L. Assessment of therapeutic efficacy of miR-126 with contrast-enhanced ultrasound in preeclampsia rats. Placenta. 2014 Jan;35(1):23-9.
13. Lu CL, Wang Y, Yuan L, Li Y, Li XY. The angiotensin-converting enzyme 2/angiotensin (1-7)/Mas axis protects the function of pancreatic β cells by improving the function of islet microvascular endothelial cells. Int J Mol Med. 2014 Nov;34(5):1293-300.
14. Ding H, Hong C, Wang Y, Liu J, Zhang N, Shen C, Wei W, Zheng F. Calreticulin promotes angiogenesis via activating nitric oxide signalling pathway in rheumatoid arthritis. Clin Exp Immunol. 2014 Nov;178(2):236-44.
15. Tang Y, Yin Y, Miao L, Wei B, Zhai K, Ji G. Nitric oxide enhances

- extracellular ATP induced Ca^{2+} oscillation in HeLa cells. *Arch Biochem Biophys.* 2015 Jan 1;565:68-75.
- 16. Cao J, Xie H, Sun Y, Zhu J, Ying M, Qiao S, Shao Q, Wu H, Wang C. Sevoflurane post-conditioning reduces rat myocardial ischemia reperfusion injury through an increase in NOS and a decrease in phosphorylated NHE1 levels. *Int J Mol Med.* 2015 Dec;36(6):1529-37.
 - 17. An XZ, Zhao ZG, Luo YX, Zhang R, Tang XQ, Hao D, Zhao X, Lv X, Liu D. Netrin-1 suppresses the MEK/ERK pathway and ITGB4 in pancreatic cancer. *Oncotarget.* 2016 Apr 26;7(17):24719-33.
 - 18. Wu KW, Kou ZW, Mo JL, Deng XX, Sun FY. Neurovascular coupling protects neurons against hypoxic injury via inhibition of potassium currents by generation of nitric oxide in direct neuron and endothelium cocultures. *Neuroscience.* 2016 Oct 15;334:275-282.
 - 19. Yin S, Zhang S, Tong G, Deng L, Liang T, Zhang J. In vitro vasorelaxation mechanisms of Isoapiole extracted from Lemonfragrant Angelica Root on ratthoracic aorta. *J Ethnopharmacol.* 2016 Jul 21;188:229-33.
 - 20. Wu Q, Zhong ZM, Zhu SY, Liao CR, Pan Y, Zeng JH, Zheng S, Ding RT, Lin QS, Ye Q, Ye WB, Li W, Chen JT. Advanced oxidation protein products induce chondrocyte apoptosis via receptor for advanced glycationend products-mediated, redox-dependent intrinsic apoptosis pathway. *Apoptosis.* 2016 Jan;21(1):36-50.
 - 21. An XZ, Zhao ZG, Luo YX, Zhang R, Tang XQ, Hao D, Zhao X, Lv X, Liu D. Netrin-1 suppresses the MEK/ERK pathway and ITGB4 in pancreatic cancer. *Oncotarget.* 2016 Apr 26;7(17):24719-33.
 - 22. Wu KW, Kou ZW, Mo JL, Deng XX, Sun FY. Neurovascular coupling protects neurons against hypoxic injury via inhibition of potassium currents by generation of nitric oxide in direct neuron and endothelium cocultures. *Neuroscience.* 2016 Oct 15;334:275-282.
 - 23. Yin S, Zhang S, Tong G, Deng L, Liang T, Zhang J. In vitro vasorelaxation mechanisms of Isoapiole extracted from Lemonfragrant Angelica Root on ratthoracic aorta. *J Ethnopharmacol.* 2016 Jul 21;188:229-33.
 - 24. Wu Q, Zhong ZM, Zhu SY, Liao CR, Pan Y, Zeng JH, Zheng S, Ding RT, Lin QS, Ye Q, Ye WB, Li W, Chen JT. Advanced oxidation protein products induce chondrocyte apoptosis via receptor for advanced glycationend products-mediated, redox-dependent intrinsic apoptosis pathway. *Apoptosis.* 2016 Jan;21(1):36-50.

Version 2016.12.08