



碧云天网站



微信公众号

碧云天生物技术/Beyotime Biotechnology

订货热线: 400-168-3301或800-8283301

订货e-mail: order@beyotime.com

技术咨询: info@beyotime.com

网址: http://www.beyotime.com

## Propidium Iodide/碘化丙啶

产品编号	产品名称	包装
ST511	Propidium Iodide/碘化丙啶	5mg

### 产品简介:

- Propidium Iodide简称PI, 中文名为碘化丙啶。分子式为 $C_{27}H_{34}I_2N_4$ , 分子量为668.40, 纯度>95%。
- 进口分装, 常用于细胞凋亡(apoptosis)或细胞坏死(necrosis)的检测, 常用于流式细胞仪分析。

### 包装清单:

产品编号	产品名称	包装
ST511	Propidium Iodide/碘化丙啶	5mg
—	说明书	1份

### 保存条件:

4°C避光保存。

### 注意事项:

- 本产品对人体有刺激性, 操作时请小心, 并注意适当防护以避免直接接触人体或吸入体内。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

### 使用本产品的文献:

1. Wu G, Liu ZS, Qian Q, Jiang CQ. Effects of Berberine on the Growth of Hepatocellular Carcinoma Cell lines. Medical Journal of Wuhan University. 2008 Jan;29(1).
2. Liu Y, Sheng Z, Liu H, Wen D, He Q, Wang S, Shao W, Jiang RJ, An S, Sun Y, Bendena WG, Wang J, Gilbert LI, Wilson TG, Song Q, Li S. Juvenile hormone counteracts the bHLH-PAS transcription factors MET and GCE to prevent caspase-dependent programmed cell death in Drosophila. Development. 2009 Jun;136(12):2015-25.
3. Li DL, Liu JJ, Liu BH, Hu H, Sun L, Miao Y, Xu HF, Yu XJ, Ma X, Ren J, Zang WJ. Acetylcholine inhibits hypoxia-induced tumor necrosis factor- $\alpha$  production via regulation of MAPKs phosphorylation in cardiomyocytes. J Cell Physiol. 2011 Apr;226(4):1052-9.
4. Cao X, Deng W, Wei Y, Su W, Yang Y, Wei Y, Yu J, Xu X. Encapsulation of plasmid DNA in calcium phosphate nanoparticles: stem cell uptake and genetransfer efficiency. Int J Nanomedicine. 2011;6:3335-49.
5. Meng LY, Liu HR, Shen Y, Yu YQ, Tao X. Cochinchina momordica seed extract induces G2/M arrest and apoptosis in human breast cancer MDA-MB-231 cells by modulating the PI3K/Akt pathway. Asian Pac J Cancer Prev. 2011;12(12):3483-8.
6. Zhao Q, Xue Y, Wang JF, Li H, Long TT, Li Z, Wang YM, Dong P, Xue CH. In vitro and in vivo anti-tumour activities of echinoside A and ds-echinoside A from Pearsonothuriagraeffei. J Sci Food Agric. 2012 Mar 15;92(4):965-74.
7. Tu Z, Ma Y, Tian J, Li H, Akers W, Achilefu S, Gu Y. Estrogen receptor  $\beta$  potentiates the antiproliferative effect of raloxifene and affects the cell migration and invasion in HCT-116 colon cancer cells. J Cancer Res Clin Oncol. 2012 Jul;138(7):1091-103.
8. Zhou Z, Wan Y, Zhang Y, Wang Z, Jia R, Fan Y, Nie H, Ying S, Huang P, Wang F. Follicular development and expression of nuclear respiratory factor-1 and peroxisome proliferator-activated receptor  $\gamma$  coactivator-1 alpha in ovaries of fetal and neonatal doelings. J Anim Sci. 2012 Nov;90(11):3752-61.
9. Jun Fang, Meihu Ma, Yongguo Jin, Ning Qiu, Chan Wang, Guodong Ren and Xin Huang. Assessment of Salmonella enteritidis Viability in Egg White during Early Incubation Stages by Fluorescent Staining Method. Asian Journal of Animal and Veterinary Advances. 7: 556-67.
10. Zhen-Jun S, Yuan-Yuan Z, Ying-Ying F, Shao-Ju J, Jiao Y, Xiao-Wei Z, Jian C, Yao X, Li-Ming Z.  $\beta$ ,  $\beta$ -Dimethylacrylshikonin exerts antitumor activity via Notch-1 signaling pathway in vitro and in vivo. Biochem Pharmacol. 2012 Aug 15;84(4):507-12.
11. Tu Z, Li H, Ma Y, Tang B, Tian J, Akers W, Achilefu S, Gu Y. The enhanced antiproliferative response to combined treatment of trichostatin A with raloxifene in MCF-7 breast cancer cells and its relevance to estrogen receptor  $\beta$  expression. Mol Cell Biochem. 2012 Jul;366(1-2):111-22.
12. Feng C, Xu Z, Li Z, Zhang D, Liu Q, Lu L. Down-regulation of Wnt10a by RNA interference inhibits proliferation and promotes apoptosis in mouse embryonic palatal mesenchymal cells through Wnt/ $\beta$ -catenin signaling pathway. J Physiol Biochem. 2013 Dec;69(4):855-63.
13. Hou Y, Chu M, Du FF, Lei JY, Chen Y, Zhu RY, Gong XH, Ma X, Jin J. Recombinant disintegrin domain of ADAM15 inhibits the proliferation and migration of Bel-7402 cells. Biochem Biophys Res Commun. 2013 Jun 14;435(4):640-5.
14. Li Q, Zhou X, Shi Y, Li J, Zheng L, Cui L, Zhang J, Wang L, Han Z, Han Y, Fan D. In vivo tracking and comparison of the therapeutic effects of MSCs and HSCs for liver injury. PLoS One. 2013 Apr 30;8(4):e62363.
15. Zhou R, Huang W, Yao Y, Wang Y, Li Z, Shao B, Zhong J, Tang M, Liang S, Zhao X, Tong A, Yang J. CA II, a potential biomarker by proteomic analysis, exerts significant inhibitory effect on the growth of colorectal cancer cells. Int J Oncol. 2013 Aug;43(2):611-21.
16. Wang Y, Jiang XL, Peng SW, Guo XY, Shang GG, Chen JC, Wu Q, Chen GQ. Induced apoptosis of osteoblasts proliferating on polyhydroxyalkanoates. Biomaterials. 2013 May;34(15):3737-46.
17. Yang F, Huang W, Li Y, Liu S, Jin M, Wang Y, Jia L, Gao Z. Anti-tumor effects in mice induced by survivin-targeted siRNA delivered through polysaccharide nanoparticles. Biomaterials. 2013 Jul;34(22):5689-99.
18. Zhou S, Wu H, Zeng C, Xiong X, Tang S, Tang Z, Sun X. Apolipoprotein E protects astrocytes from hypoxia and glutamate-induced apoptosis. FEBS Lett. 2013 Jan 16;587(2):254-8.
19. Nie C, Yang D, Liu N, Dong D, Xu J, Zhang J. Thyrotropin-releasing hormone and its analogs accelerate wound healing. J Surg Res. 2014 Jun 15;189(2):359-65.