

## DHA (≥98%, Reagent grade)


产品编号	产品名称	包装
ST1267-250mg	DHA (≥98%, Reagent grade)	250mg
ST1267-1g	DHA (≥98%, Reagent grade)	1g

### 产品简介:

CAS Number	Chemical Formula	Molecular Weight	Purity	Grade
6217-54-5	C <sub>22</sub> H <sub>32</sub> O <sub>2</sub>	328.49	≥98%	Reagent grade

### ➤ 基本信息(General Information):

Name (Chinese)	DHA
Name (English)	DHA cis-4,7,10,13,16,19-Docosahexaenoic acid
Specifications	Reagent grade, ≥98%
Chemical Formula	C <sub>22</sub> H <sub>32</sub> O <sub>2</sub>
Synonym (Chinese)	顺式-4,7,10,13,16,19-二十二碳六烯酸
Synonym (English)	DHA
Beilstein Registry No.	1715505
EINECS Number	-
MDL Number	MFCD00065722
UNSPSC Code	12352211



### ➤ 产品描述(Description):

Biochem/physiol Actions	二十二碳六烯酸DHA是一种具有22个碳原子和六个双键的ω-3多不饱和脂肪酸，第一个双键出现在甲基末端的第3位(22: 6 n-3)。DHA是脂质膜和髓鞘的组分。DHA还可作为前列腺素和类二十烷酸等信号分子的前体。
-------------------------	--

### ➤ 性质(Properties):

assay	≥98%
storage temp	-20°C
shipped in	dry ice

### ➤ 安全信息(Safety Information):

Hazard Pictogram Codes	-
Signal Word	-
Hazard Statements	-
Precautionary Statements	-
Personal Protective Equipment	Eyeshields, Gloves, half-mask respirator (US), multi-purpose combination respirator cartridge (US)
Hazard Codes (Europe)	-
Risk Codes (Europe)	-
Safety Codes (Europe)	-
RIDADR	NA 1993 / PGIII
WGK Germany	3
RTECS	-
Flash Point (F)	143.6 F
Flash Point (C)	62°C

### 包装清单:

产品编号	产品名称	包装
ST1267-250mg	DHA (≥98%, Reagent grade)	250mg
ST1267-1g	DHA (≥98%, Reagent grade)	1g

—	说明书	1份
---	-----	----

**保存条件：**

-20°C保存，五年有效。

**注意事项：**

- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

Version 20190710