

吡啶菁绿(Purity \geq 95%)

产品编号	产品名称	包装
Y026024-100mg	吡啶菁绿(Purity \geq 95%)	100mg
Y026024-500mg	吡啶菁绿(Purity \geq 95%)	500mg
Y026024-1g	吡啶菁绿(Purity \geq 95%)	1g
Y026024-5g	吡啶菁绿(Purity \geq 95%)	5g

产品简介:

CAS Number	Chemical Formula	Molecular Weight	Specifications
3599-32-4	C ₄₃ H ₄₇ N ₂ NaO ₆ S ₂	774.96	Purity \geq 95%

➤ 基本信息(General Information):

Name (Chinese)	吡啶菁绿
Name (English)	Indocyanine Green
Specifications	Purity \geq 95%
Chemical Formula	C ₄₃ H ₄₇ N ₂ NaO ₆ S ₂

➤ 性质(Properties):

Storage and Sensitivity	RT
Appearance	-

包装清单:

产品编号	产品名称	包装
Y026024-100mg	吡啶菁绿(Purity \geq 95%)	100mg
Y026024-500mg	吡啶菁绿(Purity \geq 95%)	500mg
Y026024-1g	吡啶菁绿(Purity \geq 95%)	1g
Y026024-5g	吡啶菁绿(Purity \geq 95%)	5g
-	说明书	1份

保存条件:

室温保存, 三年有效。

注意事项:

- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。